

## Product datasheet for **SC215147**

### **SLC37A3 (NM\_032295) Human 3' UTR Clone**

#### **Product data:**

Product Type:	3' UTR Clones
Product Name:	SLC37A3 (NM_032295) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	SLC37A3
ACCN:	NM_032295
Insert Size:	1541 bp



[View online »](#)

**Insert Sequence:** >SC215147 3'UTR clone of NM\_032295  
 The sequence shown below is from the reference sequence of NM\_032295. The complete sequence of this clone may contain minor differences, such as SNPs.  
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTCAACACTGCCAGCCACCTGAGATGCTGACCAGCAGTGAAGGCTGCACTAGTTTTCTACTACAGGA
ATTACTGTTGATGATTTTTGCTGTTGTTTCATGAATGTGCCGAGTGACTTCTGACTCTGCCAGCATCTT
TGACCTCTGACTCTGCCCGCTCTTTGTTTGGCCTATTGAGGTTTAGCGTTTTAACCTTAACTGTTGTT
TTGGACTTGAACCTTCTCTAGCCTTAAAAATGGAGTGAACAGGAGGCAAAGCCACGCGTGGTTTGTGTT
GGTGGTATGTTCCACAAATGGAAGAGATCATTGCATTTTGTCTTGGTAATGAGTTTACAGAGAGTGGGT
CATCTTTGAAGCATCAAATCCTGATCAAATGCTCAAGCTGGACCGAGCTGCAAGCCTGCCTTGGACAAA
GAATCCCTCTGTTCACTGTGGTGGGAAAGATCACACTACATGCCTGTTGATTGGCTCAGTCACTCTGTG
TCTGATCTAAACGTCATTAGCCCTAGAAGCATGAATTGCTTCAAATTATTGTCAACTTGTCTCTTCC
ATTTTCATTCAAATTAACCTTTGACTCCTGGAAGTATTGAGTCTTCTTTCAAGGACCATAAGGTACATC
AGTTATCTTGAACATTCTGACGTGTAAAGGACCATAAGGTACATCATTTATCTTGAACATTCCGACGTG
TAACTCACAGCCGAAGCACTCCATCCCAGATTTGTTGTCTGGGAACCTTAGGATCCTCTAGGAAGCTAAT
CTGCTTAGTCTATTTTTAGAGGATTGATCTCTGGCACAAGCAGCTTCTGGAGTTACCTCAGCAGAAGA
CTAGAATTAGAGAAAAGGAAAGACCTTTTCTTAGTAAGTGTCAAGTACAGATGCTCTTTGACTTAT
GATGGGGTTATGTCCAGTAAACCCATTGGAAGTCAAATGCATTTAATACCCCTGACTGGACACCACAG
CTTAGCCTCGTCCACCCTAAATGTACTCAGAACACTTAGAATTGCCACAGTTGGGCAGAATCATCTAA
CACAAACTCTATTTTATAATCGAGTGTGTAATATCATGTAATTTATTGAATATTGTACATTATGTTG
AAATTGCAACCATTTTACACCATTGTAAAGTCCAAAAATCAGAAGCTGAGGACTGTCTGTACATGACAT
TGAATCACATGAATGGATTACTAGGGCCTGAGGGGTGTCGGTTCTTTAGAAAAATGTAGAGTATATATA
TCTTCAGGGCTGTGAATCCATACTAGTTATCAATATTTTAAAGTCACTTGTATATACATATTTTAAA
TTATGCAAGCAAAGTAGGAGGAGAATGTGCAATTTTGGAGGATACCTAATTTGCATTGCGTTAGGGGAT
ATTTTTCAACCTCTTGTCTTATACTCTGATATGGGGTGTGTTTATGAACTTGTTAATGTTATTTAATCA
TGCAGATTTCCAATCAGTGTATTGCTCATTACTGACAAAAAGGGAAAAATACGAAACACAGAAAGGGAA
AATAAACTGACTTTTATAATGAA
AGCGGACCGACTTACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCC
CAACCTGCCATCACGAGATTTGATTCCACCGCCGC
  
```

**Restriction Sites:** Sgfl-RsrII

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_032295.4](#)

**Locus ID:** 84255

**MW:** 59.8