

Product datasheet for **SC215168**

ZNF780A (NM_001142578) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF780A (NM_001142578) Human 3' UTR Clone
Symbol:	ZNF780A
Synonyms:	ZNF780
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001142578
Insert Size:	1525 bp



[View online »](#)

Insert Sequence: >SC215168 3' UTR clone of NM_001142578
The sequence shown below is from the reference sequence of NM_001142578. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site
Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**GCGATCGC**

ACACAGGTGAGAAGGCATCT**TGA**ATGTAAGGAATGTGGAAAGTCCTTTAATCGTGTCTCAAACCATGTTCAACATCAGAGTATTCGTGCTGGTGTAAAACCAT**TGA**ATGTAAGGGGTGTGGAAAGGCTTCATTGTGGTTCAAACGTTATTCAGCATCAGAAAATTCATTCCAGTGAGAACTCTTTGTATGCAAGGAGTGGAGGACGACCTTTAGATATCATTACCACTTATTTAACATTACCAAATTCACACTGGTAAAAAACCTTTGAATGTAAAGAATGTGAAAGGCCTTCAGTTTTCTGACACAGCTTGCTGGACTTCAGGTCAATCATACTGGCGAGAAGCCATTTGAATTTGAGAAATGTGGGAAGGCCTTCAATAGTGGCTCAAACCTTGTTCAACATCAGAGTATTCATACTGGTGAGAACTCTA**TGA**ATGTAAGGAGTGTGGGAAACTTTCAGATTTAGTTTGGTCCTGACTGCACATCAAAGAATTCATACTGGTATGAAACCTGTGAATGTACAGAATGTGGGAAGACCTTTAGTTGTAGCTCAAACATTGTTCAACATGTTAAATTCATACTGGTGAGAAACGCTATAATGTAAGGAATATGGGAAAGCATTTGCTGTGGATGATCAGCTGTCTCAACATCAGAAAATTCAGAATAGTTAGAAAATTTGTGACAATAAGGAGTGTGGACAAGCCTTCACTGTGTACGAAAAATTTACTCAACACCAGAGAGCTCATATTGATGAGAAATCTTGTAACATAAAAAGATGTGTGAAGAACTTTTCATCATGGCCTGGGATTTGCTCAACTTCTGAGTATTCATCTTCTGAGAAATCTTTGAGTGAAGGAATGTGGGGCGTTTTATAGCCACACTAAATACCTTAGAGTTCCGAGAAGATAATTCCTGTGAGAAAGTCTTTGAATGAGGAATATAGGGTAGTTGTCATGTTCCAAATGTATCCCTATGAAATAATTTACCAGATACAGGCATACCTCGCTGAATTGTGATTGCAGTTATTGTGTTTTATACAACTGAAGGTTTTGTGGCAACCCTGTGGTGAGCAAGACGATCAGGTCCATTTCCCAACGTCAGGTGTGATTTTCATGTCTTTATGTCACATTTTGGTAATTCTCAGAATATTTCAAACATTTTCATTGTTATATCTATCATGGTGATCAATGACCAGTGATTTTGTGTTGCTATTGTAATTGCTTTGGGATACCACAAATCTGCCCCATAAAGATGGTAAACCTAGTTGATACTTGTGTGGGTTTTGAAAGAATAGTCAAATATACAATCATATACCCCAACCTAGATTCAAATTTGTTAACAATTTGCACGTTATTTTTGTATTTGTGTACCATTGAGAATAAGTTTTAGACATCATGACACTTTACCCTTTATTTAATCAACGTGTGTATCAACCAACATGAGACAATATAGACAATATCATATAGCCATTGTACCATAATTACCGTTAAGGAAAG

ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCG

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: NM_001142578.1

Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Locus ID: 284323