

Product datasheet for **SC215170**

ZNF780A (NM_001010880) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZNF780A
Synonyms:	ZNF780
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_001010880
Insert Size:	1525 bp



Insert Sequence:

>SC215170 3' UTR clone of NM_001010880

The sequence shown below is from the reference sequence of NM_001010880. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site **Blue**=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**GCGATCGC**

ACACAGGTGAGAAGGCATCT**TGA**ATGTAAGGAATGTGGAAAGTCCTTTAATCGTGTCTCAAACCATGTTCAACATCAGAGTATTCGTGCTGGTGTAAAACCAT**TGA**ATGTAAGGGGTGTGGAAAGGCTTCATTGTGGTTCAAACGTTATTCAGCATCAGAAAATTCATTCCAGTGAGAAACTCTTTGTATGCAAGGAGTGGAGGACGACCTTTAGATATCATTACCACTTATTTAACATTACCAATTCACACTGGTAAAAACCTTTGAATGTAAAGAATGTGGAAAGGCCTTCAGTTTTCTGACACAGCTTGTGGACTTCAGGTCATTCACTAGCGAGAAGCCATTTGAATTTGAGAAATGTGGAAAGGCCTTCAATAGTGGCTCAAACCTTGTTCACATCAGAGTATTCATACTGGTGAGAACTCTA**TGA**ATGTAAGGAGTGTGGAAACTTTAGATTTAGTTTGGTCTGACTGCACATCAAAGAATTCATACTGGTATGAAACCTGTGAATGTACAGAAATGTGGGAAGACCTTTAGTTGTAGCTCAAACATTGTTCAACATGTTAAAAATTCATACTGGTGAGAAACGCTATAATGTAAGGAATATGGGAAAGCATTTGCTGTGGATGATCAGCTGTCTCAACATCAGAAAATTCAGAATAGTTAGAAATTTGTGACAATAAGGAGTGTGGACAAGCCTTCACTGTGTACGAAAAATTTACTCAACACCAGAGAGCTCATATTGATGAGAAATCTTGTAACATAAAAGAATGTGTGAAGAACTTTATCATGGCCTGGGATTGCTCAACTTCTGAGTATTCATACTTCTGAGAAATCTTTGAGTGTAAAGGAATGTGGGCGTTTTATAGCCACACTAAATACCTTAGAGTTCGGAGAAGATAATTCCTGTGAGAAAGTCTTTGAATGAGGAATATAGGGTAGTTGTATGTTTCAATGTAATCTCCCTATGAAATAATTTACCAGATACAGGCATACCTCGCTGAATTGTGCATTGCAGTTATTGTGTTTTATACAACTGAAGGTTTGTGGCAACCTGTGGTGAGCAAGACGATCAGGTCCATTTTCCCAACGTCAGGTGTGATTTTCATGTCTTTATGTACATTTTGGTAATCTCAGAATATTTCAAACATTTTCATTGTTATATATCTATCATGGTGATCAATGACCAGTGATTTTTGATGTTGCTATTGTAATTGCTTTGGGATACCACAAATCTGCCCCATAAAAGATGGTAAACCTAGTTGATACTGTTGTGGGTTTTGAAAGAATAGTCAAATATACAATCATATACCCCAACCTAGATTCAAATTTGTTAACAATTTGCACGTTTATTTTTGTATATTGTTGTACCATTGAGAAATAAGTTTTAGACATCATGACACTTTACCCTTTATTTAATCAACGTGTGTATCAACCAACATGAGACAATATAGACAATATCATATAGCCATTGTACCATAATTACCGTTAAGGAAAG

ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCG

Restriction Sites:

SgfI-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NM_001010880.2](#)

Summary:

May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Locus ID:

284323