

Product datasheet for **SC215171**

SKT (KIAA1217) (NM_001098501) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SKT (KIAA1217) (NM_001098501) Human 3' UTR Clone
Symbol:	SKT
Synonyms:	RP11-324E23.1; SKT
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001098501
Insert Size:	1568 bp



[View online »](#)

Insert Sequence:	>SC215171 3'UTR clone of NM_001098501 The sequence shown below is from the reference sequence of NM_001098501. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC ACCCCATCCACAGCAAAAGAAACCTCTTAAGGTCAAATCCTATTAGGCACAAGTCGGAGTTACATTTA AAAAAATTAACAGTCTACAACAAGTGTTCACAAGAGAATGTAACATATTGCTGTATCGTTTGAGGC TTAATGCTAAATATGTGCTAAATACTGGATTAATAGATTTAGTAAAGCTCGTTCGTTTGTGTTT TCTTTTACCTAGTTGCTATAGTGTCTACAGTCTATACTCAATACCTATAAAATGCAGTAAGCATGTGT TACAGAAAGAGGTTCTGGTGGGAGAGAAAGGTGCGTGTGAGACAGGAGAATTGCTTAAGCATATAAAA CATGTATGATTCCAGAATTTAGTATGTTTGTATAAACTATTTTTCATTACGGAGACTAGAAGTGAA CAGAGAATTACACAAGTGTGACTATACAAATTGTAACAGATACTATAATTTCTTTTATTTTATG GTTATTTAGCTTTATTACAGATTCTATTTTGTCAAACTTCATGGTTCCTTTCAAGATCTTTTGTG CAAAACATTTTGATACTATAGCATTGTACATTTGAAAGTAGTGTCTAGACTATAAAACCAATGAACCT CTACATGAGCCCTACAGACAGGCATGTGTAGAAGGCAATTTATCAAACCTATTGCACTGCCATGAAAAG TGTGTATAATAATTTGCTAGCCCAAGCAAGCTAGTTTCTTTGCTTGCTTCTTTCTTTCTTTTCTT TCCTTTTTTTTTTTTTTTTTCTTTTTAACATGTTGAGATTCTCTAGTTGTTTCTTTGGCGTATCTAAC CCCTTCTTTTGTCTGAGACCTGGTAACCCACGCTCTTGCAATTGTGGATTTTAAATGTATACTCTG TACGGTTCTGTAAACCGAAAACTTTTGTAAATATATAATATACATAGACATAAAAACTGTATGTG ACAGCACATAGAGTAGTTTCCACACCAAGTAAATTTTATGCATGCTTTAAAGTATATATCGGGA CCGGCAGAAATGGAAGTATCCATACATTTTAAAAAGCAACAAGTTGCACAGCTAGAGTGTGTTTGT AATAAATGTATTTGTATAACACAGTCATGTAATATACAGAACTATAAGCAGAGACTTTGCAAACTAAA TAAAGGGCTGCATGCTTATTATTTTGTACCTTGCTACTATACTACTCTAGTCAAAGAACGAAAT GTAAGTGTACCGAGTTAAATGTTTTCCGCTTGAGGGATGTAACCACATCCACTCAGAGGACACTAC TTTTCTGAAAGCTCTGGGTGACTAATGATGAGTTCCTAATAAATTAATTGCAAGTGTGGTGCCTTGGA TGTGGCTGTGGCTCGCTTCTCTGTGGCTATCAAGGTGTAGATGACAGAAAGCAAACTGGAT ACAGAGTTTCCACCTCAGTTCCTGGAGGGCTCTATTATTTTCTCTTTTTAAAAAAGTCCAGTA GAAGTAAAGTGGAAATAAATGTCTTTATCAAAAAAAAAAAAAAAAAAAAA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001098501.1</u>
Summary:	Required for normal development of intervertebral disks.[UniProtKB/Swiss-Prot Function]
Locus ID:	56243
MW:	61.3