

Product datasheet for **SC215175**

FAM190A (CCSER1) (NM_207491) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	FAM190A
Synonyms:	FAM190A
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_207491
Insert Size:	1569 bp



Insert Sequence:

>SC215175 3'UTR clone of NM_207491

The sequence shown below is from the reference sequence of NM_207491. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
TTCAAGACCACATACCCTAGGGAATTCGAAAACCTGACTAACATCAATTGAAGATTGGGGTTCAGGAG
GCCTAGATGATCCACCATGACTTATTCAAAAGTTGAAGATGAAAGGAAAGAATAAAGACTCTTTCCATC
CAGTTCTCTAGTCTGAAAGGAAAAGAAATAGGCAAAGGACAGTAGGCTTTTCTATATCATTATATTCAT
TCAAATGATTCTCTTAAATAATAATGAACTAAGGAGGCTCTGTAAGCAAGACAAACGGTCAGGAGGCC
TCCTTGCTCCAAGAATACTCCCTTTTACTTTTATTAATGTTTTGTAAGGTAAACTAAATTGACTAA
AGTATCAATATTTCCATATATGCGTGATCTTTAGAGACCATAAACTGTCAGTGTATTAGGCTGTCTT
GCATGGCTATAAAGAAATACCTGAGACTAGGTGATTTATAAGAAAAGGGGTATAATTGGCCACAGTTC
TGCAGACTGTACAGGAAGCATAGTAGCATGTGCTTCTGGGAGGCCTCAGGAAGCTTCCAATCGTGGCG
GAAGGCAAAGGGAAGCAGGCACATTGCATGGTGAGAACAGAAACAAGAAAGGAGGTGGCTGGGAG
GCGCCACATAATTTTAAATTACCAGATCTCACAAGAACTCACTATCCGGAAGACAGCACCCAGCCATGA
GGGATCCACCCTCATGATTCAAACACCTCCCACCAAGCCCCACCTCCAGCACTGGGGATCACAACCTCGA
CATGAGATTTGAGCAGAGACAAATTTCCAACTATATCACTCAGAACTGGAATAGAAAGAGTATGCAC
CATAACACTTTTATTCTCTACAGAAACAGACACTGTGCTGAGCTCAAAAAATAGTGATTCATAAATAT
GAATGTGCAAACCTCTTTTGTACTGTATCTGTGGAATCCTGGTTACTGGGGCTTTGAGCCCCGTGCCCA
TAAATCTCTGTAATACTAACCTAAGTATGATGATTTCTTTAAATCAAATTTAATGAAATGATGGACA
TAAAGCCCATGGAAGTGTGTAATCTCTCTTATTTGATTAGATGTTAACTTGAACACTTTTTT
TTGATAAAATATGAAATACATAAATGAGTTAAATGTTGTGACTTTCCCTATAAGATTGAAAACAATAA
GGTGTTTTTAATACTTTTTTCTGTGTTCTCTTGAAAGTAGCTTTACATGTATTAAGCAATCCTGGCTT
TAGATATTAGGACAATAAAAAATCTGCAATTTCCAGTTTTATATAGCCATTGTTTGCAAGTATGA
AAACCTGCTTTTTGTTACTGAAATATTATCTTTATTTGGCTTTTTGATACTTTTATGTTTGCCTT
GAGCAAGTCACTTATCTTTCTATTTGTGTTTTCTATTTTCAGTTTATTTATGAAAGACAGTAATAATGA
TTGTGATGTACCTTACAAAAAAGATAAAATCGTATTGAAGAGACTGTTGAATATGAACGACGTGTATT
AAATGCCCGATATAGGAACAAAAATAAAATGAAGAACATTTCTATTATAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:

SgfI-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

NM_207491.2

Locus ID:

401145

MW:

61.3