

Product datasheet for **SC215279**

SH3D20 (ARHGAP27) (NM_001159330) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SH3D20 (ARHGAP27) (NM_001159330) Human 3' UTR Clone
Symbol:	SH3D20
Synonyms:	CAMGAP1; FLJ35572; FLJ43547; MGC120624
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001159330
Insert Size:	1588 bp



[View online »](#)

Insert Sequence:	>SC215279 3'UTR clone of NM_001159330 The sequence shown below is from the reference sequence of NM_001159330. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC CAGTGCGCGGACATCTTCCCGCCGACACTGACTGCTGGCCTGTGACTGGGGCGGCGCCGCGGTCTCTGCC ACACAAGCTGGGCGGCGGAGGCCACGCAGCCGGCCTTCTTCTCTCTGGGACCCCTCCGCCAGCGCATAG CCGAGGCCGGTGTGACTTCTGCACCCTCGGTTCTGAGGGTACGGTGACCCCTAGTGGGCAGTTTGCAA AATGTGATTCTTCTTCCCAACTCCCCATCCCCCTTCCCTTCCCGTCACGTCTGTTTGGGGTTAAT TCGGTTTTTCTCTGTTGCATCGCGCTACTGTGCGTGTGCGATAGCGTGTGTTGGGGTGAGAGTTTGT TTTCTGGAATGGTAGGTGCTGGGAGGAGAGTTTGTGAGGGGCTTCTGGCTGCTTCTGGCCCTCACC TCGTGGAGGCCCTTACAGAGACCCTGTGGGCCCTGGCCCTGTGCTGGCACTGTGCCAGTCATGAGGCAG CTCTGATCACTTCCCCACTGTGGAACAGGACTGACCCAGCCTTCAGTGTGGGCTGCTGAAGCTATCCT CCTCAGGCCTCAGGGATGACCTCCTGCCTGAGCCTCTCACAGGCTGGCTGTGGGCCAGTTTCATCTGCT TTCCTGTTGGGGTCCCAGGCCCTGCTGTCCTTGACCACTGGTGTCTGTGCAAGGCTTCTTCCCAT TCACCAAGTGACACCTTGCATCTGCCGCTCGGCATGCACCAAGTCCACACACCATCCCATTTTACAGA CAAGGACGCTGAGGCCTGCAGCAGCAGTGTGACTTGCTCAAGGTCCAGTGAGTGACCTCATTCCCAGA AAAGGCTCCTCCACACCAGAGTACAGCCTGGGTAGGGGAAAATCAGTTCTTTAGCTACCAACCATC CAACCTTTGGGCCTATGTGAAAAGAAAGAACTAAGCTGGGTGTGTTCTGTCTGGACCTGGGGAGGCC CTGAAGGCAAGAGGGAACTGTCCAGCTGTTCTGTCTAGGGGAGGGGACATAGCCCTAGCAGGAG CTCCAGCCCCCTTGGCACTGTGACACACAAGTACACCATCTGGGGCCGCTTTGCCACGAAGAGCT GGGCAGGCCTGCAGGGTGTGGGGAAGGAGGACACAACCTCAAGAAAGGAAGCGTGAACCCAGGGAACA GCGGGTCCCTTCCCTCCTCAGACACAAGCCACCTCAGCTTGTGGCTTGTGGCCCCAGCCCCACCAACC CACCTGTTCAATTATTCAACAGACAATGACAGCTGATATTTATTGGACATTTGCACCATGCCAAGCATT CGGCTTGATTATCCCATTTGTTTCTCACAGCCGGTATTTATTGTCTGCTCCTCTGTGCCAGGTGCTGT GCTCTGGGCAGGGGCACTGCATGGGCTGCCTGCCCTGGTGGAGCTTGTGGTCTGATGGGTGAGGCTGAC CCAAGCCCAACCCATTGCCAACAGGGCCAGGGCAAGAGTACACACAGGGGCTCATACCATATGTCTAA ATATTTAAAGTTATCAATCAAGCTAACAAGTGTAAATAAAATATGTTCTATTCTCTACTTTGAAAAA ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001159330.1</u>
Summary:	This gene encodes a member of a large family of proteins that activate Rho-type guanosine triphosphate (GTP) metabolizing enzymes. The encoded protein may pay a role in clathrin-mediated endocytosis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Aug 2013]
Locus ID:	201176

MW: 56.4