

Product datasheet for **SC215281**

LPHN3 (ADGRL3) (NM_015236) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LPHN3 (ADGRL3) (NM_015236) Human 3' UTR Clone
Symbol:	LPHN3
Synonyms:	CIRL3; CL3; LEC3; LPHN3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_015236
Insert Size:	2000 bp



[View online »](#)

Insert Sequence: >SC215281 3'UTR clone of NM_015236
The sequence shown below is from the reference sequence of NM_015236. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

[illegible]

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_015236.7](#)

Summary:	This gene encodes a member of the latrophilin subfamily of G-protein coupled receptors (GPCR). Latrophilins may function in both cell adhesion and signal transduction. In experiments with non-human species, endogenous proteolytic cleavage within a cysteine-rich GPS (G-protein-coupled-receptor proteolysis site) domain resulted in two subunits (a large extracellular N-terminal cell adhesion subunit and a subunit with substantial similarity to the secretin/calcitonin family of GPCRs) being non-covalently bound at the cell membrane. [provided by RefSeq, Jul 2008]
Locus ID:	23284
MW:	77.2