

Product datasheet for **SC215314**

Myelin oligodendrocyte glycoprotein (MOG) (NM_206814) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: Myelin oligodendrocyte glycoprotein (MOG) (NM_206814) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: MOG
Synonyms: BTN6; BTNL11; MOGIG2; NRCLP7
ACCN: NM_206814
Insert Size: 1210 bp
Insert Sequence: >SC215314 3'UTR clone of NM_206814
The sequence shown below is from the reference sequence of NM_206814. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
TTCCTTGAAGAGCTACGAAATCCCTTCTGAGTGATGTACATCTTGGCAGGGGTGGAGGAGAGCCTGGT
TGCCCAGGGATTTGTCTTGGGGACATCTCATCCATCAAGTTGCACACTACTGGCATCTTTGCTATGG
GGACATTTCCAATTTGCACTTTCAGGAACACTCTGAATTTCCAAGTAGAATTGATTTCCCTTCTTGTCA
TCTACCTTTTCTTCTTATTTCCATTTTTATTACCCTTCTTTCCATTTCTCTCCAGTCTTCCACCT
GGAAGCCCTCTCTGGCTAAGGACAGGCAGGTGCCCTCTCTCCATCAGAGGACACCTGTAAGGAGAGC
AACACAGGATGGTCTCTGCCATGAAGTGGAGGCCAGGAATCTCCTCACTGAAAATTACAGTATGGTAAC
TTTGCAAATGGTGGTGGTTTCTTCCAAGACTCCAGCCCTGATTGCGCAAAACTGAAAGGCATGTGAAGG
GAAGGAAGAGGAAGAGTGAACAACTTGAAGAGAGAGCTGAGTGAGCTGAAGAGTGAGGATATGAGTAG
CCCCAACCCAAACCTGGAGATGGGGAGAAACCTACAGAATACTAGCCAGAGCTCCTCCTTGTCTTGCA
GCCTACTAGGGACCTGGGAAGCAAAAACGAAAGCTGGGCAACATGCCTGCTTTAGAATGTTTTCCCTTC
TACTTACACATCTTCCACAGGTCTCAGAATCTTTCCTTCTCTCATCCTTTTCTCTATCTTCATATCT
ATCAGAGTATCCACTGTTTATTTCAACAATACTACTTGTGATGGTGCAGACACAAACAAGCTAGGTGC
TAATTAATAAAGATACGAGTTTTGGCCGGTGCGGTGGCTCACGCCTGTAATCCCAGCACTTTGGGAGG
CCGAGGCGGGCGAATCACGAGGTCAGGAGTTCAAGACCAGCCTGGCCAACATGGTGAAACCCCATCTCT
ACTAAAAATACAAACAATTAAGTGAATAGTGGTGGGCACCTATAATACCAGCTACTCCGGAGGCTGA
GGCAGGAGAATCGCTTGAACCCAGGAGGCGAGAGTTGCAGTGAGCTGAGATCGTGCCACTGCACTCTAG
CCGGAGTGACAGAGTAAGACTCTGTCTCAAAAATAAATAAATAAATAAATAAATAAATAAATAAATAA
TAAATAAAAAATAAATAACAAGTTTTTCATAAGCACA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_206814.6
Summary:	<p>The product of this gene is a membrane protein expressed on the oligodendrocyte cell surface and the outermost surface of myelin sheaths. Due to this localization, it is a primary target antigen involved in immune-mediated demyelination. This protein may be involved in completion and maintenance of the myelin sheath and in cell-cell communication. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]</p>
Locus ID:	4340
MW:	44.9