

Product datasheet for **SC215322**

NIBAN2 (NM_001035534) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	NIBAN2 (NM_001035534) Human 3' UTR Clone
Symbol:	NIBAN2
Synonyms:	bA356B19.6; C9orf88; FAM129B; MEG-3; MINERVA; OC58
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001035534
Insert Size:	1532 bp



[View online »](#)

Insert Sequence: >SC215322 3'UTR clone of NM_001035534
The sequence shown below is from the reference sequence of NM_001035534. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

GGCAAGTTGGACGCCCGCAAGATCCGCAGATTCTCATTAAAGGCCAAGGAGGGCGGAAAGATCGCGGT
TAACAATTGGCAGAGCTCAGAATTCAAACCGATCGCC
GACAGTGCAGGGGTGCAGACTGAGTTCTAGCGCAGTGGGTCCCTGACTGCTGCACATGGCACAGGCCGT
TCCCTTCCGGACCCAGGCAGGCTCAGCTCTGGGGAGGGCACCTGGTCTGTGCCTTGTGGGTGGAGCG
GGGCAGGGCTGTGTGCCACCCGAGGGAGCGGGGCCACCTGAGTCACTTTATTGGGTTTCACTCAACAC
TTCTTGTCTCCCTGTTTTCTTCTGTGGGATGATCTCAGATGCAGGGGCTGGTTTTGGGGTTTTCTCTGC
TTGTGCCAAGGGCTGGACACTGCTGGGGGGCTGAAAAGCCCTCCCTTCTGTCTTCTGTGGCCTCCA
TCCCCTCATGGGTGCTGCCATCCTTCTGGAGAGAGGGAGGTGAAAGCTGGTGTGAGCCAGTGGGTTT
CCGCCCACTACCCAGGAGCTGGCTGGGCCAGGACCGGGAGAGGGAGCACTGCTGCCCTCCTGGCCCTG
CTCCTTCCGCAGTTAGGGGTGGACCGAGCCTCGCTTTCCCACTGTTCTGGAGGGAAGGGGAAGGAGGG
GGTCTTCAGGCTGGAGCCAGGCTGGGGTGCTGGGTGGAGAGATGAGATTTAGGGGTGCTCATGGGG
TGGGACGGCCTGGGGTGAAATTGAAAAGGCCAGACGTCGAGGTCGCGAGGGGAAGTGTCTGAGT
GAAGGAGGGGACCCCATCCTGGGGGATGCTGGGAGTGAGTGAGTGAGATGGCTGAGTGAGGGTTATGG
GGAGCCTGAGTTTTATGGGCTGTGTATCCCTTCTCCCGGCCCGAGCCTGCTCCCTCCTGCCCGCC
TGGCCACAGGTCTCCCTCTGGTCCCTGTCCCTCTGGTGGTTGGGGATGGAGCGGCAGCAAGGGGTGTA
ATGGGGCTGGGTTCTGTCTTCTACAGGCCACCCGAGGTCCTCAGTGGTTGCCTGGGGAGCCGACGGG
GCTCTGAGGGGTACAGGTTGGTGGGCCCTCCCTGAGGGTCTGGGGTCAAGCTTTGGCCTCTGCTGCC
TCTCAGTCACCAAGTCACCTCCCTCTGAAAATCCAGTCCCTTCTTTGGATGTCCTTGTGAGTCACTCTG
GGCCTGGCTGTGTCCTCCTCAGCTTCTTGTTCCTGGGACAAGGTCAGGCCAGGATGGGCCAGGCC
TGGGATCCCCACCCAGGACCCCCAGGCCCTCCCTGCTGCTTTGCGGGGGGAGGGCAGAAATGG
ACTCCTTTTGGGTCCCCGAGGTGGGGTCCCCTCCCAGCCCTGCATCCTCCGTGCCCTAGACCTGCTCCC
CAGAGGAGGGGCCCTTGACCCACAGGACGTGTGGTGGCGCCTGGCACTCAGGGACCCCCAGCTGCCGAG
CCCTGGTCTCTGCCGATCTCTTCCCTCTTGTCCCGAAGATCTGCGCCTCTAGTGCTTTTGGGGGT
CCCATCATCTCTCCCTGATATTGTATTGAAAATATTATGCACACTGTTTATGCTTCTACTAATCAATAA
ACGCTTATTTAAA
ACCGGTGAGCGCGCGGCATCTAGATTGAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCGCCTTCTATGAAAGG

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: NM_001035534.3

Summary: May play a role in apoptosis suppression. May promote melanoma cell invasion in vitro.
[UniProtKB/Swiss-Prot Function]

Locus ID: 64855

MW: 53.4