

## Product datasheet for **SC215327**

### **CREM (NM\_182853) Human 3' UTR Clone**

#### **Product data:**

Product Type:	3' UTR Clones
Product Name:	CREM (NM_182853) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CREM
Synonyms:	CREM-2; hCREM-2; ICER
ACCN:	NM_182853
Insert Size:	1595 bp



[View online »](#)

**Insert Sequence:** >SC215327 3'UTR clone of NM\_182853  
The sequence shown below is from the reference sequence of NM\_182853. The complete sequence of this clone may contain minor differences, such as SNPs.  
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAAGATCTTTATTGCCATAAAGTAGAGTAACTGTCTTTGACTTGGACCTTGTTTACTCTAATCAAGGCA
GGAGATGCAGCAGTCTACTTATTGCCATGTGGACTTGTGGGAAGGACACGTGTGACCCTTAAGAATCC
AGTTTGGATTAGTGTGTTGAAATTGAATTGGGAATGTTGTTCCAGGATGTGGAATGCAGCCGTGATCACA
CTTACCGAGCTTACTTTGATCTGTTTGTCAATAGCATGCAAAAAATGCTTTGTTGGCCCTTTGCTTCTG
CTTTTTTTCAGGGAAGCTGCCAAAGAATGTCGACGTCGAAAGAAAGAATATGTAAATGTCTGGAGAGC
CGAGTTGCAGTGTGGAAGTCCAGAACAAGAAGCTTATAGAGGAAGCTTAAAACCTTGAAAGACATTTGT
TCTCCCAAAACAGATTACTAGAAATATTTAACTATGAACTGAAGGCAGCATGTATAGTTGCTTTTGAAG
GAATACAATATATAGCTGGCAAGATGGTGGCTTCTTTTCTTTGTATCATTTCATCTTCTTCTTAATCA
CTTAACATTCTAAAATGCTTCACTGTACGTAGTTAAGTCGTAGCTATAAAGCTTCAAATTTTTTAAAGA
GACAAACTGTAAAAATGTGTGATTCTTAAAATGCAATATTTGTAAGGCTTGTCCAATGCCACATAC
TTGCAGCTCCCATTCTATGTGTCATCAATAGTGTCTATGCAATAAAAATTTTGACGGTCTTTAAATC
ATTTTAGGAAAGGATGATCAAAAATAATGCATCCAGCAGTACAATAAAGTAAACCACAAAAAATACC
TCAGGAAAGAATAGAAAGAAAGTCTATCTAATGACATGCCTATATGAGAAGAATAGCCTAGATATGAAT
ATATGGCATTGTCAGATTTTATATTAGTTGCTTTGTTAAAAAAGATTGTATTGCTGTCCTTGA
ATGCCATAGTCAAAGAGAGTTTTAATAGAACCATGTTGGTTGCACCTTTGTAGTGTGTTGGTGCCTATT
AAATATCTGAACATTTACTACAGTTTTTAACTCTACTGTGTAATATAAAGATCTTGCAGAAGTTCTTA
GTGTTGGTTTAAATTACCTAATGAAAGTATGACATATTTTTTATATCTGGAATGAGCCTGTTGGGAT
CGCATTGCATACCTTCGGGTACAAAGTCAAGTCTTACTACTGGGTGCTGATCCTTGACGCGCCCTTCT
ACCATCTCACGGTGGGATGGCCGACGGCTGTGCACCCAGAAGAAGATGGCTGTAGCGTTTCGGCCT
TCATAATGGCCTGAGACTTTCTTCTGTAGGTGGTCTGGAGCTGTCCGGCTGGTGGCCCCCTATTTTGC
CATTTAGCGAACAACCACAGGAATTTAAAAACAAAAACATCCCAAGATTTTTTCATTTCAAATGCTT
CAAAGTCCACATTAGATCAGATACTCCGCTGTCCGCACATTCAGCTGAGGTTTCATTACAATCGAGACTG
CAATGTGATCTATGTTTCATCTGTTTTTATAATAAAAAGCTTCAGGGAACAAGCCCAAAGCCCTCAC
ACAAAAA
ACGCGTAAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** SgfI-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_182853.2](#)

**Summary:**

This gene encodes a bZIP transcription factor that binds to the cAMP responsive element found in many viral and cellular promoters. It is an important component of cAMP-mediated signal transduction during the spermatogenetic cycle, as well as other complex processes. Alternative promoter and translation initiation site usage allows this gene to exert spatial and temporal specificity to cAMP responsiveness. Multiple alternatively spliced transcript variants encoding several different isoforms have been found for this gene, with some of them functioning as activators and some as repressors of transcription. [provided by RefSeq, Jul 2008]

**Locus ID:**

1390

**MW:**

62.1