

Product datasheet for **SC215357**

TMEM130 (NM_152913) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TMEM130 (NM_152913) Human 3' UTR Clone
Symbol:	TMEM130
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152913
Insert Size:	1596 bp



[View online »](#)

Insert Sequence: >SC215357 3'UTR clone of NM_152913
The sequence shown below is from the reference sequence of NM_152913. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TATAAGTCTGTCAAAACTTACACCGTGTAGCACTCCCCCTCCCCACCCCATCTCAGTGTTAACTGACT
GCTGACTTGGAGTTTCCAGCAGGGTGGTGTGCACCACTGACCAGGAGGGTTTCATTTGCGTGGGGTGT
TGGCCTGGATCATCCATCCATCTGTACAGTTCAGCCACTGCCACAAGCCCTCCCTCTGTACCCCT
GACCCAGCCATTACCCATCTGTACAGTCCAGCCACTGACATAAGCCCACTCGTTACCACCCCTT
GACCCCTACCTTTGAAGAGGCTTCGTGCAGGACTTTGATGCTTGGGGTGTTCGTGTGACTCCCAGG
TGGGCTGGCTGCCCACTGCCATTCTCTCATATTGGCACATCTGCTGTCCATTGGGGTTCTCAGTT
TCCTCCCCAGACAGCCCTACCTGTGCCAGAGAGCTAGAAAGAAGGTCATAAAGGGTTAAAAATCCATA
ACTAAAGTTGTACACATAGATGGGCACACTCACAGAGAGAAGTGTGCATGTACACACACCACACACAC
ACACACACACACACACACACACAGAAATATAAACACATGCGTCACATGGGCATTTTCAGATGAT
CAGCTCTGTATCTGGTTAAGTCGGTTGCTGGGATGCACCTGCACTAGAGCTGAAAGGAAATTTGACCT
CCAAGCAGCCCTGACAGGTTCTGGGCGCGGCCCTCCCTTTGTGCTTTGTCTCTGCAGTTCTTGCGCC
TTTATAAGGCCATCCTAGTCCCTGCTGGCTGGCAGGGGGCTGGATGGGGGCGAGGACTAATACTGAGTG
ATTGCAGAGTGCTTTATAAATATCACCTATTTTATCGAAACCCATCTGTGAACTTTCACTGAGGAAA
AGGCCTTGACGCGGTAGAAGAGGTTGAGTCAAGGCCGGGCGCGGTGGCTCACGCCTGTAATCCAGCAC
TTTGGGAGGCCGAGGCGGGTGGATCACGAGATCAGGAGATCGAGACCACCTGGCTAACACGGTGAAC
CCGCTCTCTACTAAAAAATACAAAAAGTTAGCCGGGCGTGGTGGTGGGTGCCTGTAGTCCAGCTACT
CGGGAGGCTGAGGCAGGAGAATGGTGCAGAACCCGGGAGGCGGAGCTTGCACTGAGCCAGATGGCGCCA
CTGCACCTCCAGCTGAGTGACAGAGCGAGACTCTGTCTCAAAAAAAAAAAAAAGAAGAGGTTGAGTCAG
CAGGGACTTGGGTTCCCTGTGTGTGAGGGGGCATTCTTGCTGCCAGCTGCTCCGAGGTGGCCTTGA
GAAGGAAGAAGCAGGATGACAGAGCCTGAGCAGCGGAACAGCCTGCACCCTCCCTCTGCCCCAGCGA
CCTGGGCTGTGGCTGAGACAATAATGAGGCCAGAAGTAGCCGGAGCCTGTGAGGAAGGGCAGGGGAGGA
CTGTGGGGTCTGGGCTCTGTCGCTGTAACCATCTGCTCCAGGCTGTGTGCAGAAAATGGCATTTACAC
TATTGTGCAGCTATTCTCATGAAATACTGCCATTGTTGCTAAATAAAGCTTGTGTGCTCTGAATATAG
CCAAGCATA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_152913.3](#)

Locus ID: 222865

MW: 59.6