

Product datasheet for **SC215358**

TMEM130 (NM_001134450) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TMEM130 (NM_001134450) Human 3' UTR Clone
Symbol:	TMEM130
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001134450
Insert Size:	1596 bp



[View online »](#)

Insert Sequence:	>SC215358 3'UTR clone of NM_001134450 The sequence shown below is from the reference sequence of NM_001134450. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TATAAGTCTGTCAAAACTTACACCGTGTAGCACTCCCCCTCCCCACCCCATCTCAGTGTTAACTGACT GCTGACTTGGAGTTTCCAGCAGGGTGGTGTGCACCACTGACCAGGAGGGGTTTCATTTGCGTGGGGGTGT TGGCCTGGATCATCCATCCATCTGTACAGTTCAGCCACTGCCACAAGCCCTCCCTCTCTGTCAACCCCT GACCCAGCCATTACCCATCTGTACAGTCCAGCCACTGACATAAGCCCACTCGGTTACCAACCCCTT GACCCCTACCTTTGAAGAGGCTTCGTGCAGGACTTTGATGCTTGGGGTGTTCGTGTGACTCCCAGG TGGGCTGGCTGCCCACTGCCATTCTCTCATATTGGCACATCTGCTGTCCATTGGGGGTCTCAGTT TCCTCCCCAGACAGCCCTACCTGTGCCAGAGAGCTAGAAAGAAGGTCATAAAGGGTTAAAAATCCATA ACTAAAGTTGTACATAGATGGGCACACTCACAGAGAGAAGTGTGCATGTACACACACCACACACAC ACACACACACACACACACACACAGAAATATAACACATGCGTCACATGGGCATTTTCAGATGAT CAGCTCTGTATCTGGTTAAGTCGGTTGCTGGGATGCACCTGCACTAGAGCTGAAAGGAAATTTGACCT CCAAGCAGCCCTGACAGGTTCTGGGCGCGGCCCTCCCTTTGTGCTTTGTCTCTGCAGTTCTTGCGCCC TTTATAAGGCCATCCTAGTCCCTGCTGGCTGGCAGGGGGCTGGATGGGGGCGAGGACTAATACTGAGTG ATTGCAGAGTGCTTTATAAATATCACCTATTTTATCGAAACCCATCTGTGAACTTTCACTGAGGAAA AGGCCTTGACGCGGTAGAAGAGGTTGAGTCAAGGCCGGGCGCGGTGGCTCACGCCTGTAATCCCAGCAC TTTGGGAGGCCGAGGCGGGTGGATCACGAGATCAGGAGATCGAGACCACCTGGCTAACACGGTGAAC CCCGTCTCTACTAAAAAATACAAAAAGTTAGCCGGGCGTGGTGGTGGGTGCCTGTAGTCCCAGCTACT CGGGAGGCTGAGGCAGGAGAATGGTGCGAACCCGGGAGGCGGAGCTTGCAGTGAGCCAGATGGCGCCA CTGCACCTCCAGCTGAGTGACAGAGCGAGACTCTGTCTCAAAAAAAAAAAAAAGAAGAGGTTGAGTCAG CAGGGACTTGGGTTCCCTGTGTGTGAGGGGGGCTTCTTGCTGCCAGCTGCTCCGAGGTGGCCTTGA GAAGGAAGAAGCAGGATGACAGAGCCTGAGCAGCGGAACAGCCTGCACCCTCCCTTCTGCCCCAGCGA CCTGGGCTGTGGCTGAGACAATAATGAGGCCAGAAGTAGCCGGAGCCTGTGAGGAAGGGCAGGGGAGGA CTGTGGGGTCTGGGCTCTGTGCTGTAAACCATCTGCTCCAGGCTGTGTGCAGAAATGGCATTTACAC TATTGTGAGCTCATTCTCATGAAATACTGCCATTGTTGCTAAATAAAGCTTGTGTGCTCTGAATATAG CCAAGCATA ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001134450.2</u>
Locus ID:	222865
MW:	59.6