

Product datasheet for **SC215409**

SIRPB2 (NM_001122962) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SIRPB2 (NM_001122962) Human 3' UTR Clone
Symbol:	SIRPB2
Synonyms:	dj776F14.2; PTPN1L; PTPNS1L3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001122962
Insert Size:	1619 bp



[View online »](#)

Insert Sequence:	>SC215409 3'UTR clone of NM_001122962 The sequence shown below is from the reference sequence of NM_001122962. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC ACCTTAGCATGGAGCAAGGGTCAAGAGTAGGGGTCAGCCCCAGAGTGAGGACCCTCTGAGTTGGAGAG GAGCCAGGGCTCCTCAACCATTTCCTACCTCCAGTCCCAGCCTCTAGGTGCCCCAGGCCTCATGACA AACTCCTAGATCCCTACATCTGGTTTTGGTCCACCTAGTGAAATTCCTTCTTTCACCGGGCTTCCT CTAAAATGTCTCCCTTTCTTTTTGGCCTGTTCAAGACCTCCTTGCTTTTCAGTCCCTGGCTCAGTCT CTCCTCAACACCCTTGCCCTGCTGCAGCCCTTTCTGGTGCGCCCTGCCCTTTCCACCTCGTACA TCCTTCTTGGCTCCAACATCCAACCTCAGAGTCTTCTCCAGGAGATGTCTGTAAGAATCTCTGAAC CAACCAGCCAGACCATCTGTGCCCTCCATCTACACCTTTCTCCCACTCCTTCTGCCTTCCTTCCAT CCCCCTCATGGCTGGCTTGGGCAGGTATAATATTAGAATGCAGGTTGAGCAACTATAACAAAGCTCTTA AATAACAGTGGCTTAAACAGTGGAAATCAACCAGAAAGTTGACCATCAGCAGGCCAAGCAATACAGAG ACTCCCTGGTATTGAGACCCAGGATTCACTGATCTCATTGCTACCAGGTCCACCTTCTAGGCAGCCAGA CTGGAAAAGAGGGCAGGAAAGGGGAGCAGGACCCTCCCCTTTAAGTGCACAGTCAGGAACCTTGGCCACC TCATTTATCTCTACTTGGCTGGAATGTGGTCACATGGTCACACCTAGCTGCAAGAAACACTGGGAGATG TAGTCTTTATTTCTGGCAGCAATGCGCCAGCTGCAAGTTTTCTACTAGAGAAACAGATGGCAGATATC AGGGGATAACCAAGTTATCTCCACCACAGCAGCATACAGACAGCCTCTCACCTGCCCTGTGGGACACCTG AGTTCAATGCCAGCTAGCTAGCCAGCACTTCTTCCCACTATCACCTCCCCTGGGGCAGCATGATGTGG GGCAGTAGTTCCCAAGATGAGTGATTTTGCCCCACTGGACTTTTGGCAATGTCTAGAGATGTTTTGG TTGGCACAACTGGGGGGGTGCTACCACCATCTAGTGGACTGAGAAGCCCTGACATGGGGAAGAGTGTG CATGCCAGGAGTCAGACACACCTGCCTTTAACCTGAGGCCTCTGCCTCCTCCTGTGCACCTCAGT GACTAATCAGAGTCCCTTCCCATCACGGAACATCCAGGATACTAATGTGGACTTCTCTGCATTGTGTAA GAACCAATTCAAGACCAGGCACGGTGGCTTATGCATGTAATCCCAGCACTTTGGGAGGCCGAGGTGGGT GGATCACCTGAGTTCAGGAGTTTGGAGCCAGCCTGGCTAACATGGTGAAACCTCGTCTCTACTAAAAAT ACAAAAAATTAGCCAGGCGTGGTGGTGTGCACCTGTAATCCCAGCTACTTGGGAGGATGGGGCAGGAGA ACCGCTTGAACCTGGGAGGCAGAGGCTGCAGTGAGCTGAGATCGCGCCATTGCACTCCAGCCTGGGCAAC AAGAGCAAACTCCGTCTCAAAAAAAAAAAAA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001122962.2</u>
Locus ID:	284759
MW:	58.9