

Product datasheet for **SC215421**

PTPN20B (PTPN20) (NM_001042360) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PTPN20B (PTPN20) (NM_001042360) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	PTPN20
Synonyms:	bA42B19.1; bA142I17.1; CT126; PTPN20A; PTPN20B
ACCN:	NM_001042360
Insert Size:	1611 bp



[View online »](#)

Insert Sequence: >SC215421 3'UTR clone of NM_001042360
 The sequence shown below is from the reference sequence of NM_001042360. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GTTCTTCGGAAACTTCTGACTTTGGATTAAAGAAAGACTTCTGTTGCCTCTCACTTGAATTACCAAGTG
GGTTTGCACCTCCTCATAAAGAACATGTTTGCAGTGTGCTGAAGGGCTTTGCTATGCATACAATCTGCT
TTCTTGGTTTATCAGTTTATTTCTTTCTAAAAGCTCCCTGAAGGGCAATATCATTTGGCTTGGGGTGA
TCAGTGTTTACTTATTGATCTTGTAGACAATATCAAAATAAATTCCCACATTTTCCAGTGAAACAGAT
GTTACATAAAAACGATTGCAGCTTGGCTATTTGGTTGAAGGGATTACAGAGCCCAATAAAGGATTTAAAA
TATATTCATTAAGATTTTATTTGGAAAGGTGGCTGGAGAGAGCTGAGGATTTCCAGGACTTTGTAAGTT
CTTATTCTGGGAGAACATAAGGCCAATAATCATGACCTCTCCAGGCATTTTAAAGACAGATGTCTATT
CATGTTCTTTAGCTAGAGCCTGTACTTTTTGCTGGCATTGTAATAACCCAGTTTAAAAAGAGTCCAGTT
AGGGTGGACTAACTTTGGACACAATTGGCTTCCATTTCTACATTTTCACTACTGCTGCCTTCTACAG
CTGCTAGACCAAGACCTGTTGGTCTGGGAAGCATTTCATGGATAGGGAGAGCTCCTCTCGGTGAACAGT
CCAAAATAAAATAGATGTTTATATAGAAAAGCCCAAGAGGAGACTTTTGCATGCCTGAGTTCTTTCT
ATCCCACCTAACACTTAACATATTACTTAGTCTGCTTTGTTAAAAGCAAGTATTACCTTAACTTGCC
TCTTACTCTTTGCCCTTTAGCTAACTAATAAAGTTTGATATAGGCATTATTATATAAATTCTGAGTCATT
CATGGTATCTCTCATGTTTGATGATTTTTCAAATAAGATCTATGATAGTTTTTTTTCCAGAGTTCCA
TTAAATCATTTATTTCTTTACTTTCTCACCTCTGTTGAAACATTTAGAACTGGATTTGGGAACCCAA
TTTTGGAAAACAGATTCATAGTCATGAAAATGGAACTTCCATATTCTGTTTTTGGAAAGATGTGGCC
ATTATTACAGTAATTTTATTATAGGACTTTGCCTCGTACAATTAATAGTGATATTTTGGACAAGGAGTT
CTGGTGACAAGCTATACCTAATTATAAGCTATAAAACAATAGATATGAGTGTGTTGTACAGTTAACTCA
ATGGAGATCAGAAATTTCTATGTATTGAGAAAATGTTAATATCAATCTATAAACTTTGAATTTCTAAG
AGGCTTATTTTGTCTTTGGCTGAATGAGTATATTTGAATTGGTTGAATAATTAATAATTCTCATTGT
AAAAATAATTATATGCCAAAATATATTTGATGTTAAATCAAATAGATGATTCTGTTTACATTGTTTCA
ATGAATAAATAATCTGTGTTAATTTCAATTTTATAATTGGCCTTTAATATTGTATCTCTAATTTTATTT
TCTCTCTGTTACTGTAATAATAAGCTATAATGTATAACAATTTTCTTCAGAAGAATTCTATGCTATTA
TTAAAATAAAATATTTACTGTGTA
ACGCGTAAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001042360.4](#)

Summary: The product of this gene belongs to the family of classical tyrosine-specific protein tyrosine phosphatases. Many protein tyrosine phosphatases have been shown to regulate fundamental cellular processes. The encoded protein appears to be targeted to sites of actin polymerization. A pseudogene of this gene has been defined on chromosome 10. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014]

Locus ID: 26095

MW: 62.8