

Product datasheet for **SC215428**

CCDC43 (NM_001099225) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CCDC43 (NM_001099225) Human 3' UTR Clone
Symbol:	CCDC43
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001099225
Insert Size:	1597 bp



[View online »](#)

Insert Sequence:	>SC215428 3'UTR clone of NM_001099225 The sequence shown below is from the reference sequence of NM_001099225. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC CCGAAACACCAATGTGGAAGATGTCCTTAATGCCCGAAAACCTGGAGCGAGACTCACTTCGGGATGAATC CCAAAGGAAGAAGGAACAGGACAAGCTGCAGAGGGAGAGAGACAACTAGCCAAGCAGGAGCGCAAGGA AAAGGAAAAAAGGACACAGAGAGGGGAGCGAAAGCGATAACCTTGGCCTACTCTATTCTTTCCCCC CATGAACCTATCAAAGGGAGAACTGGTTGTGTATTTAAGAGAAATGAGAGTGCATTGCAAAGTAGTTTC CCAAGGCATTTCCCTATTGTTTACATGGTCAAAAGGAAAATCAAATCACTTCTAATTACAAAATGTG CTGTTTTGGTGGGTGGCAATCAGATTATTATAGTTGATGACTGTACCAAGATCTGGAATGGGGGCC ACCTTTATATTTAATACTCAAGTCCTGTGCTCCTTTTGGCCACTTTAAAAAGTCTTCCCTCACCTTGG ATAAATAGGATTAGAAAGCTAACAGTTTTTCATTCTGAAGGAGCAGAAACATAATCAATCCCTCCCTGA GATAAAATGTGGGCATGCACATTGCTTTGGAGTCAGTATTCAGCTGTAAGTGATGATTGTTTTATGGC TGAAATAGAAAGCTAATGAAGAGGCCTACCTGGATCCAGATTTTATGTTAACATGGCAGCAAACTAGG AAAATGGGAAGAAATTACTTGGGGGTGTTCTAGATCTGATTTTTTTTTTTTTTTTGAAGTTTCAGCC TGTGGGAGGAAAAGAGGTATCTTTAAATAGAGATAGTTACATATTTTCTCATCTAGGCTGTGTAC TATTGGGTGCAAGGATATCATTGTTGGCTAGATTCAATTTGTATAATCATGTATCCTCTTGTGTGCTGG TAGAGATTTTAATCCTGATTTTCCATAAAACATGAGTATTAAGAAATAATCCTGGTTTGGAGAACT GGAGAAATCACCTTTTAAGGAAGAAACACTGGAATTTCTGCTAACACCAAGATATTTAAGAGTGTAC ATAGTAGGTGCTCAACAAATTTATTGAATGAATGAGTGAATGGAAAACTGGGAGAGTCAAAAGTGAGC AGAAGCTCTCCATTTCTACTTCTGTACAAACCACATTAATTTGTAATAAGGCCCTTCTCCACTTGAC TTCAGGCAGCAGATTGTCTAGAAGCCTAAGGACAGCAATTTCTCTGACAAGACAAAGTAGATATTTTAT ACCAGGGGTTGGCAAACTACTGCCACGGGCCGAATTTGGCCAGTCTGTTTTGTATGGTGCAAACTA AAAATGATTTTTACATTTTAAAGAGTTATAAAAGAAAAAATATGTGGTCTGTGAAATCTAAAATATT TACTACCTGGCCTGTTGGAGGAAAGTTTGCCAATCTCTGTTTTATACCATTAACATGAGATTACAA AACTTTTACCTTTGTGCAGAAGTTAAAAAATCAATGGTTAAGGAAAAGGAGACATGTTACCT CTTCATACACTCCTATACTGTGGCATTGCAAAAAATAAAATAACACCTTTAAAAAATAAATCTTA TTAAATTGC ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001099225.2</u>
Locus ID:	124808
MW:	61.7