

Product datasheet for **SC215519**

CTL5 (SLC44A5) (NM_001130058) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CTL5 (SLC44A5) (NM_001130058) Human 3' UTR Clone
Symbol:	CTL5
Synonyms:	CTL5
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001130058
Insert Size:	1628 bp



[View online »](#)

Insert Sequence:	>SC215519 3'UTR clone of NM_001130058 The sequence shown below is from the reference sequence of NM_001130058. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC GAGGAAAATCCACAACTAGGAAGCAGTAGAAGACAACTGGTCGTCCTACAGCTGTGTGTACCTTT TCTCCATCTGCTGTGTCTGTGCAACATTTGTTTCATAAGTGCTTTGTGTTTAGCAACACTGTATTACG ACCTTGTGGCTTGCATTTGCATGTTTTATACCAAGCTTATACTGTACTATGTGAAGCCATCAGAAGT CGCAAGGGAATTGTTAATAACATAAAACATTTTATACTAAGATCATTTGTTTGTAAATCGTTTTTAA AGAGTGGCTTGGATGTTTGAATAACTACTGAATATGTTAATATTCTTTAAATCTTAGATTGAAAAA TGATACATTACTTAAATTGATAGCTCCTAATATATTTTTAAAATTACAATAAAAGAAGACTTCTCTG CAGGGAATTTGGTCAGCAAGTGAAATTAATAATTTTAAAGTTTTCCCACTCTCGTTGGACAGTAAA TCAGTGAAAGGACTGCCCGAGTTGAGAGTTTGCTCTCTTTAAGTATAGAATGTTTCCTCTTAAACAAAT TGCCAATCATCCAGCCTTTACTACTTAGCCCTCTGACAAAGTGCCCTACTGGCTATTTAATATTACCA GCTTTTATGGGCAAGTTTACAAACATTGTTTTTAAAAAATTAACCTGCAATGTTTCGTGATTAAAA CAAGTCTTATTGCATTTGTTTCACTCTTAGCTCACTGATTGAAAAACATTTGTCATTTTGTCTGTTTG ATATCCTCACTATTATGGAATACATTGTGCAGCTAAACAATTTCCCTTGCCTAGTGGACATTCATGA ATGTGTACTACACGCAAGAAGAAACAAACCCCGAAAGAACACTTGTGGATTTCCTTTGTTTTTTTTT ACTAAAAGAGAAGTTTTAAATGAAATGTTTTCTATAGTAGATCTTTGAAATACAATAGGTATAATAC TGCATTTCTCAGTGTTTTACAAAGATCAGAAAGAGAACTTCTAGGAATTGCAAGGGAACTTTACTC CTCGAAAGGGTGCTCACAGATGTCTGACTGAATAGCTCCCTTTTAAATGATCATTTATTTTCATCAA AGCTGTTCTATATATGCCACTTCATTTTCTAAGTTTGGTATGAAAAATCAGTTTACTTACAGTATG TTAATTGTATTGTACTACTATAAACAGGAACATAATTTCCAATTCAGTTTTAAATAATTTTACCAGTAC TACTAAGTTTAAAGGAAATTAATTCAGTTGGTTACTCAGTGTTTGTACAGAAAGAGTCCAGAAAAGTA TTCACCCTAAGAGAATGTCAATCATATAATGATAATTTGTGAAAGCTTTGAGAATCAATCATCAGTAAG TTACTATCAGTTTATAAAATATTATCACATTTGTTTAAATGTGACTTTAGATACTTTTATGCCAAAAAT AAACTCACATGAGCACATGACAGTCTGAGCTCTATAATCAGTGTGCTTCTGCTGTGCAGAAATGTTAGA AACGTATTGTCTAAATATCTTTGATAATTAATAATGTTTAAATATTTAATGAAATTTGTTGTTACTTGT TAAATCTTTTTCTTTTAAATAAGATTTAAATAAGAAATAA ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001130058.2</u>
Locus ID:	204962
MW:	64.3