

Product datasheet for **SC215534**

GPR64 (ADGRG2) (NM_001079859) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	GPR64 (ADGRG2) (NM_001079859) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	ADGRG2
Synonyms:	CBAVDX; EDDM6; GPR64; HE6; TM7LN2
ACCN:	NM_001079859
Insert Size:	1587 bp



[View online »](#)

Insert Sequence: >SC215534 3'UTR clone of NM_001079859
 The sequence shown below is from the reference sequence of NM_001079859. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GGAAGCTTACACTTTATTGAGCAAATGTGATTCCTTTCTTCTAAAATCAAAGCATGATGCTTGACAGTG
TGAAATGTCCAATTTTACCTTTTACACAATGTGAGATGTATGAAAATCAACTCATTTTATTCTCGGCAA
CATCTGGAGAAGCATAAGCTAATTAAGGGCGATGATTATTATTACAAGAAGAAACCAAGACATTACACC
ATGGTTTTTACAGATTTCTGATTTGGTTTCTTATCTTTTCTTTTATAAGAAGGTTGGTTTTAAACAATA
CACTAAGAATGACTCCTATAAAGAAAACAAAAAAGGTAGTGAACCTTCAGCTACCTTTTAAAGAGGCT
AAGTTATCTTTGATAACATCATATAAAGCAACTGTTGACTTCAGCCTGTTGGTGAGTTTGTGTCAT
GCCTTTGTTGATATAAGCTAAATCTAGTGACCCATGTGTCAAAAATCTTACTTCTACATTTTTTTGT
ATTTATTTTCTACTGTGTAATGTATTCCTTTGTAGAATCATGGTTGTTTTGTCTCACGTGATAATTCA
GAAAATCCTTGCTCGTTCGCAAATCCTAAAGCTCCTTTTGGAGATGATATAGGATGTGAAATACAGAA
ACCTCAGTGAAATCAAGAAATAATGATCCCAGCCAGACTGAGAAAATGTAAAGCAGACAGTGCCACAGTT
AGCTCATACAGTGCCCTTGGAGCAAGTTAGGAAAAGATGCCCCACTGGGCAGACACAGCCCTATGGGCT
CATGGTTTGACAAACAGAGTGGAGAGACCATATTTAGCCCACTCACCTCTTGGGTGCACGACCTGT
ACAGCCAAACACAGCATCCAATATGAATACCCATCCCCTGACCGCATCCCCAGTAGTCAGATTATAGAA
TCTGCACCAAGATGTTTAGCTTTATACCTTGGCCACAGAGGGATGAACTGTATCCAGACCATGTGT
CAGGAAAATTGTGAACGTAGATGAGGTACATACACTGCCGCTTCTCAAATCCCCAGAGCCTTTAGGAAC
AGGAGAGTAGACTAGGATTCCTTCTTAAAAAGGTACATATATGAAAAAATCATATTGCCGTTT
TTAAAAAGGCAACTGCATGGTACATTGTTGATTGTTATGACTGGTACACTCTGGCCAGCCAGAGCTAT
AATTGTTTTTAAATGTGTCTTGAAGAATGCACAGTGACAAGGGGAGTAGCTATTGGGAACAGGGAAC
GTCCTACACTGCTATTGTTGCTACATGTATCGAGCCTTGATTGCTCCTAGTTATATACAGGGTCTATCT
TGCTTCTACCTACATCTGCTTGAGCAGTGCCCAAGTACATCCTTATTAGGAACATTTCAAACCCCTT
TTAGTTAAGTCTTCTACTAAGGTTCTCTTGATATATTTCAAGTGAATGTTGGATCTCAGACTAACCAT
AGTAATAATACACATTTCTGTGAGTGCTGACTTGTCTTTGCAATATTTCTTTTCTGATTTATTTAATTT
TCTTGATTTATATGTTAAAATCAAAAATGTTAAAATCAATGAAATAAATTTGCAGTTAAGATCTTTAA
ACGCGTAAAGCGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-Mlul

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001079859.3](#)

Summary: This gene encodes a member of the G protein-coupled receptor family described as an epididymis-specific transmembrane protein. The encoded protein may be proteolytically processed as it contains a motif shown to be a protein scission motif in some members of this family (PMID: 11973329). Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]

Locus ID: 10149

MW: 61.4