

Product datasheet for **SC215562**

MLX (NM_198205) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MLX (NM_198205) Human 3' UTR Clone
Symbol:	MLX
Synonyms:	bHLHd13; MAD7; MXD7; TCFL4; TF4
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_198205
Insert Size:	1630 bp



[View online »](#)

Insert Sequence: >SC215562 3'UTR clone of NM_198205
The sequence shown below is from the reference sequence of NM_198205. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTGCACCAATTGAAAAACCAGCTTTACGACCGTTCTTGGAAACCTGGAGAACAGCCAACAAGAGGCC
CTTGAATCTCTACGTGGCCACTGAACTGCTGGGCCCGGAGACTGGACTACAACACCTCAGACTGGTCA
GCTGGTTTCTACTTGGTGTGTTGTTTTCCAGCCCCATTTATCTTCAGCGGAGCCGGTGTGTTGTT
TTGTGAAAGCTTCTGATTAATTTATTATATTGACGATAAACTCAAACCTACCCAGCCTTCCCCCACT
CCATGGAAGTCCTTGGGATGGGCGTCTGCTCTGGACACCCCAAAGAGCTCCTGCCCTCTCAGCCCTTTA
TTCAAGCCTCAGATTTCTGCTCATGATCTACATAGATTTGAAACTGTTTCTCTGTTTGGTCTCTT
GGGCAACATTTTTGGCCCAAGTTTGGGCAACATTTGGCCCAAGTTTGGGCATTTTGGCAGTAGCTGTAT
GGGAGAAAAAGAGTAAGAGGAAATATCCACAGCCATGAAGGGTAAAGGGCACCTTGTGCCTAGACT
AGGGCTGCCTGGTCAGTCCCAGGTGAGGCCAAGGGCTTTCTGGCCATCTCAGGGAGGGGCCACCAGGTT
CCTCCCCCTCACCCCATATTCCATCACCTTCCTCCTCTGCTCTGGGTGGTAAGGGAAGCCCTCCCGGTTT
CCACAGGCTATGATGCTGCATGGCAGAGGCAGGTATAACACAGCACTACATATTGGAAATTTTTATTT
TTCTAAATACCAATGCAGTTTGTACGGTTACAATTTTGAATATTAAGTACGCTCAAAATCACCCCT
TTCTGTCAAGCATATCTTGGCCTCTCCCATGTCTCAGTGTGCTGCAATTTCTCCAGGACTTGGGGT
GGGGTGAAGCGTACAAAAGATACTTAAAGGGCTCCTGGGGTACACAAGCCAGCAGGTCCTGAGTG
AAGCCGTGGGCCCTCAAATGCTCGTTTTATAGCAACCTCTCTACCCTAGTTCTCAAATTCACCTTC
TGCTTCTCAGTGTGATATCTGGCAGTTTGTACTATCCAGAGGAAATTAATATTTTTATATAAAAT
TAAATTATAATAAATATTGCCAAATGCTTTCTTTAGCATTGTTCCAAGTCTAAATGTTAACCTCAAGC
TACTGCAATTTAGACAATGAAATGGGCTGGGTCTACCCCAAGCCACCAGCCCTCATCCTCTCTACCCAG
TGCTCTGGTTTATGCTTGTCTCTGACTGCTCTGCTTAAAGGTGAAAGTAGCAGGAACAACAACAAAG
CCAACCAAAAACAAGGTAGCCAGTGCAAGACATCTCACTCTTCTGACATCCTGCAGTCCCCACCAGTCC
TGACCGTGGGCCCTCAGGGGTCTGGGAGTGTGACGTTGTAATCTTCATCCGTCTCTATCCCACTTCC
TCCTGTGAGACAGGGAGACAAGTGAATGAGATGTCACCAGGATAAGACCACAGGGAAGCAAGAAGGAA
GAGAGCTCCACTTACAAAGAACTGCTTCTGCTCTTGGGTATCCTTCAAGTATTGCATCAGACAGCTC
TGTAGCCTGACAAGAAATAAAACCACCCGTTTTTCAGATGGGCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: NM_198205.2

Summary:	The product of this gene belongs to the family of basic helix-loop-helix leucine zipper (bHLH-Zip) transcription factors. These factors form heterodimers with Mad proteins and play a role in proliferation, determination and differentiation. This gene product may act to diversify Mad family function by its restricted association with a subset of the Mad family of transcriptional repressors, namely, Mad1 and Mad4. Alternatively spliced transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]
Locus ID:	6945
MW:	60.6