

Product datasheet for **SC215576**

TPD52L2 (NM_199361) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TPD52L2 (NM_199361) Human 3' UTR Clone
Symbol:	TPD52L2
Synonyms:	D54; TPD54
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_199361
Insert Size:	1623 bp



[View online »](#)

Insert Sequence: >SC215576 3'UTR clone of NM_199361
 The sequence shown below is from the reference sequence of NM_199361. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAGCCCTGTTCGATCCCGCACCTTTCAGCCTGTGGTTGCTTACCCGCTGCAGAGCACACGCAACC
CAGCCTCAGCATCACAGCCGAGCTCTGTTGAGCGGAGCAGCCAGCCAGGCGGATGAGCAGAGCCGGC
CCTGAGGACAGTCTGCCATCCACGCGGAGATGTGGCTGCCGCTTTCATGAATTTGAAGAACACAG
GCTTGTACACAGATGTTTTACACTCACGTTTGTAGATGAAACAGATCACTGTGCTGCTCCTCCTAGGGG
TGCAGGAAGTGGACAGGGCGGAGGGTTTAAAGAATATTGAGCCAAAGCCAGGCTCCCTTTGGGAATC
ATGTTAGCCCATCAGAATGTTGAAGGATTGAAGAGTTCTAAGCATAAAAATAAGTGGCATTCTGACTT
CTTCTCCTCCTCCTCCCTGACTCACAGAAGGAATGCAATCACCCAGCAAGCTTACCTGTTACGCAA
TTTTTTATCTCAAATGCCGAACGAGAAAAGTGTCCATTTCTGAGACCCCGAAGAAAGAACTGACCC
TCAGCAGCTGCCTGATTGTTACGCGAATCTAGCTTAAACGGAAGCAAATTCATTATTTTTAAATGCAG
TGGACTTTTCAAAAAGTTTAAATTAGGCAAGCAGCTTTCAGCCTCATAGAATATTATTTCTTTGGACTC
AAGCTGAAATACAAGCCTTACATTGCCTTATGCTTTATTTCTTTCTAATTTTTATATGTATATAGATGA
GGGTTCCTTAATGGTTGTGAGCATTGTGTGGAATTTTACACCTGGCCTGCGTGGCAGCCTCTCCAGTT
GAGGTGTTTTATGTCACGCACACTCCATCCAGTGTACAAAACCTGCTTCTTCTCAACCGTGGCAGC
TCCCGTGGCTCCTATGCCCTGCCCTAAAGGGCTCTTGAGCCTCTGGGAATGGGAGGGGCCAAGAGAAG
GAAAACCTGTCTTTAGCACCTTTAAAGAAGTGTGCCCCCTTCTCAGTGTGCTTTGATGAAACGACCCAC
TGGCCCGCTCACATTTCGTGACTCCAACCTCCTGCTTGTGACTTGGGATGAAACGACCCAC
AGGTGAGGTGGAGGGTGGGGCGTGGGCATCAGCCAGGATTGCCGTTACAGTCTTTTTCTCAGGAGCTAC
AAAGATCTCTTCTGTTACTAAATAGTCGACCCAGCAGCCTCTCTCGCACACCGGGCCCTGCATGT
CAGATGGCGTGGTCTGCAGGGGAGCTCTGTGCCTTAGTGGCTCTTGGCAGGACACTGAGGGCCTGCCT
GTGGTGTGCCGCTCTGCCACTCCCGGAGGGGAAGGGCTGCTCAGCTCAAGGTGCTGTTCCGGTAG
AGCAAGTGTCTCTGACAGCCGTGTCGCCGACAGTTCAGACACCCTTGGGGATGGCACTCCACACAG
ACAGAGATGCAGGGGCCAGGGAAGCCAGCGCTCGGTGCCCTTCGTCCAGGGTAAAATCGGCCTGTGG
GGTGTGGTGAGAAGGCAGGTTGTGCGGGTGTGACCGATGTATCTTTCCCTAAAGTTATTATAATAAT
GGTAATTTGTCAATAAAGCATTCTTTGGGGAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-Mlul

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_199361.3](#)

Summary: This gene encodes a member of the tumor protein D52-like family. These proteins are characterized by an N-terminal coiled-coil motif that is used to form homo- and heteromeric complexes with other tumor protein D52-like proteins. Expression of this gene may be a marker for breast cancer and acute lymphoblastic leukemia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 12. [provided by RefSeq, Aug 2011]

Locus ID: 7165

MW: 59.8