

Product datasheet for **SC215580**

TPD52L2 (NM_199362) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TPD52L2 (NM_199362) Human 3' UTR Clone
Symbol:	TPD52L2
Synonyms:	D54; TPD54
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_199362
Insert Size:	1623 bp



[View online »](#)

Insert Sequence: >SC215580 3'UTR clone of NM_199362
The sequence shown below is from the reference sequence of NM_199362. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAGCCCTGTTCGATCCCGCACCTTTCAGCCTGTGGTTGCTTACCCGCTGCAGAGCACACGCAACC
CAGCCTCAGCATCACAGCCGAGCTCTGTTGAGCGGAGCAGCCAGCCAGGCGGATGAGCAGAGCCGGC
CCTGAGGACAGTCTGCCATCCACGCGGAGATGTGGCTGCCGCTTTCATGAATTTGAAGAACACAG
GCTTGTACACAGATGTTTTACACTCACGTTGTAGATGAAACAGATCACTGTGCTGCTCCTCCTAGGGG
TGCAGGAAGTGGACAGGGCGGAGGGTTTAAAGAATATTGAGCCAAAGCCAGGCTCCCTTTGGGAATC
ATGTTAGCCCATCAGAATGTTGAAGGATTGAAGAGTTCTAAGCATAAAAATAAGTGGCATTCTGACTT
CTTCTCCTCCTCCTCCCTGACTCACAGAAGGAATGCAATCACCCAGCAAGCTCTACCTGTTACGCAA
TTTTTTATCTCAAATGCCGAACGAGAAAAGTGTCCATTTCTGAGACCCCGAGAAAGGAAAGTACCC
TCAGCAGCTGCCTGATTGTTACGCGAATCTAGCTTAAACGGAAGCAAATTCATTATTTTTAAATGCAG
TGGACTTTTCAAAAAGTTTAAATTAGGCAAGCAGCTTTAGCCTCATAGAATATTATTTCTTTGGACTC
AAGCTGAAATACAAGCCTTACATTGCCTTATGCTTTATTTCTTTCTAATTTTTATATGTATATAGATGA
GGGTTCCTTAATGGTGTGAGCATTGTGTGGAATTTTACACCTGGCCTGCGTGGCAGCCTCTCCAGTT
GAGGTGTTTTATGTCACGCACACTCCATCCAGTGTACAAAACCTGCTTCTTCTCAACCGTGGCAGC
TCCCGTGGCTCCTATGCCCTGCCCTAAAGGGCTCTTGAGCCTCTGGGAATGGGAGGGGCAAGAGAAG
GAAAACCTGTCTTTAGCACCTTTAAAGAAGTGTGCCCCCTTCTCAGTGTGCCTTTGCATGGGCC
TGGCCCGCTCACATTGTCAGTACTCAACCTCCTGCTTGTGTACTTGGGATGAAACGACCCAC
AGGTGAGGTGGAGGGTGGGGCGTGGGCATCAGCCAGGATTGCCGTTACAGTCTTTTTCTCAGGAGTAC
AAAGATCTCTTCTGTTACTAAATAGTCGACCCAGCAGCCTCTCTCGCACACCGGGCCCTGCATGT
CAGATGGCGTGGTCTGCAGGGGAGCTCTGTGCCTTAGTGGCTCTTGGCAGGACACTGAGGGCCTGCCT
GTGGTGTGCCCGCTCTGCCACTCCCGGAGGGGAAGGGCTGCTCAGCTCAAGGTGCTCTGTTCCGGTAG
AGCAAGTGTCTCTGACAGCCGTGTCCCGGACAGTTCAGACACCCTTGGGGATGGCACTCCACACAG
ACAGAGATGCAGGGGCCAGGGAAGCCAGCGCTCGGTGCCCTTCGTCCAGGGTAAAATCGGCCTGTGG
GGTGTGGTGAGAAGGCAGGTTGTGCGGGTGTGACCGATGTATCTTTTCTTAAAGTTATTATAATAAT
GGGTAATTTGTCAATAAAGCATTCTTTGGGGAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_199362.3](#)

Summary:

This gene encodes a member of the tumor protein D52-like family. These proteins are characterized by an N-terminal coiled-coil motif that is used to form homo- and heteromeric complexes with other tumor protein D52-like proteins. Expression of this gene may be a marker for breast cancer and acute lymphoblastic leukemia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 12. [provided by RefSeq, Aug 2011]

Locus ID:

7165

MW:

59.8