

Product datasheet for **SC215604**

TRIM32 (NM_012210) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TRIM32 (NM_012210) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	TRIM32
Synonyms:	BBS11; HT2A; LGMD2H; LGMDR8; TATIP
ACCN:	NM_012210
Insert Size:	1625 bp



[View online »](#)

Insert Sequence: >SC215604 3'UTR clone of NM_012210
The sequence shown below is from the reference sequence of NM_012210. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TACCATCTGAGAAGATATTCCACCCCATAGGGGATGAGAAATTATCAGTTTCTTCTGCTCCCAAGCCAA
CTTCCCTTCCCTTAGTTCTTGGTTGTTAGTGGCACATGCAGAATAGACTCAGCCTATGTCCTGATTCCA
GCTGGGTAGTTCTAGAACTTCAGAAGCTCCATCTTTAATGTTTTTATTTGTTATGTCCCCCTCCCCGC
TTCCACCTAAATTTAGAGCTTTAAAAGATGCACTGCCAAATAGGACACACGATGGTGTAGCTGAAG
TTTGATTAGCAATTAGGCACTTCCAAGGCTTTAGTAGAGAGAGCCACTTTAGCCCTTTGTCCATGTTT
GAAATTTGCCCTTGATTAATCCTTGATTTTTCCCATTTGGCTTTGATGCCCTTGATCCATTGTTTC
CTTCTACTATAATGTGCTTCATCTGTGACACTTTCTCTTGAAGCTGATTGGATTCACTGTGCATGCT
TCAGTGGGATCTGCTCCACCTTTCAAGTACATTTAAGACATCATATCCCGTAACATTATGTCTCAGTC
TGATCGTCTTTACCAGTATGAAAGTCATTCATTTAGTGCTACCAAGGGGATACACAAGCCCTTTAGGA
AGCAGTACCTCTCGCCTGGAGGATCTGTGCCATCTTGATTGAGAATTGCAGATGTGACAGAATGGATT
GACCCTAGTTGGTTGGTATTGATGACTTCAGCCTGGAAATGCTTGCCTTTAAAGAAGCATATATGGG
TTGGAATTATGCCAAAGCATAGGAAGCTGGGAATAAGCAAACAAATGCTGATATAGTCAGCAAATTTGG
ATAGTCTCTAGGGCTCATATTTTTCACTACCTCTCTCTTCTGGCCTGTGTCTAAGGAATTGTACAA
CATAGGCCAGGGCCAACAAGTGGAGAGGTGGACACATTTTCATGTTTACTACTAAAACAACAGCAAA
ACTATTTGGTTTGTATTCTGTGTTTTCTCAAGTCAGTACATACTATTTGGTTTCAGGATTTCTTTCCA
TTTCTATCAAGCATTAAATAATTGAGAAGTCTTCTCATCTCTGGTTCATGTCTTTAATTAATAA
TACGGAATTTTGAGATGATAAGAGGGATAACCTGCGCTATGCCATTTTGCCTTCTTATCTCACTGTGT
TCTTTTCAGAGTGTAGATATCTACCCAGGTCAAATTCAGCCAAAAGCACTATTATGTAAAGATAAAT
CCAAATCTTTCTCCAAATACATGTGTACAGAATCATCCCCAAGGTTTAGAGTCACCAATACCCTATT
CTGAGCTTCTGTGTATGTTTTATGTTTATAATTACAGTTCACAAAAAGGAAGTGTCTGTTTTTAA
GATCAGAAGTGGGTGAGAGGACCTGTGATGAGGTTCTGAGGGTTTGGTGTACTAAGTAATATATCT
TACTCTAGGTGCAGGGCTGATGGCAAGCCAAAGCAACTGCCTTACTTTGATGTAAACAAAATTTGCT
AAACTGAGCTGCTCAAAATGTTTTTATGGTAGTAGTGTCTTTGGATTGTATCATAATAGCTGCCAA
AGGATAGGTAAAGAGGTCATTAATGATGTTGAAGTAA
ACGCGTAAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_012210.4](#)

Summary: The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to cytoplasmic bodies. The protein has also been localized to the nucleus, where it interacts with the activation domain of the HIV-1 Tat protein. The Tat protein activates transcription of HIV-1 genes. [provided by RefSeq, Jul 2008]

Locus ID: 22954

MW: 62.2