

Product datasheet for **SC215644**

ZFYVE27 (NM_144588) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZFYVE27 (NM_144588) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	ZFYVE27
Synonyms:	PROTRUDIN; SPG33
ACCN:	NM_144588
Insert Size:	1625 bp



[View online »](#)

Insert Sequence: >SC215644 3'UTR clone of NM_144588
The sequence shown below is from the reference sequence of NM_144588. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GCCTCGTGTAAACCAGACCTTGAGCAAGTGAAGAGAGGGTCCAACCAGGCACCCGTCCTTGGG
ACCAGCAGTAGACCCCCACTCTCCCCACCCCTGGCCACTGTGGTGTGTGCTGGCAAATGTGGCCTG
AATGCTAGGTAGGCTTCCCCTTCTTCCCTCACTCTCTCCAGCTGGATTCTGGAGCTGTTCTCCATCCAT
GAGAGTGGCTGGCAATGGCTGCTCTCAATCCCTTGAGGGAGAAGAGCCCCCTGGAGGGCTGGCATGTTT
GTCCTGCTCTGCCTGGGACTGAGCGAGTGGACTTAGGGCTGGGCAGGCAGTAGCCACCAGAGGGCAGCA
GCGAACTAGGCCAGCCCTGACTGGGGTCTGAAGATCAGGGTCAGTGTGGCTGTGCCTGGGAATCCAGA
CCTGAGGTTGGGAAAAGAGGTTTTCTCCTGCAGGGTACTGGCCAGGCCCTCAGCCTCAGAGAGCCTG
CAGAAGGGCTTGGGAGTGCCACACCCCATCTCTGCTGATTGAATGTCCCTCCAGGCACCAGGATCTCAT
CATTTCCCATCAGAGGGTGTGGCCAGGCCTAACAAAGACCATGGGTGCTTCTAGAAACAGGGTTGAAGT
TCCCAGATTCCCCTGAGAGGAGAATGTGTATAGGAGGGTTGGCTGAGTCCTTCAGCGTTAAGTGGAGGA
AAGCTTGGGGAAGCCCCAATAGCTGGACAGACCTCGGCCCTCCCCTCGAAGACACCTCAATTCACAGACT
CTCAGCCCACACAATGCCCAAGTGTCCCAGCTCCGCTGGAGCAGCTGCAGGGCACTTGGATCACAACCT
TCTGCACCCTCTGTCCAGAGTCTAGGGCAGTCTCCACTGGCCCAGCACTCCAGTTTCTTTCCCTGCC
TCTTGTCCAATGGAGTGGGAGGCCAGGTGAGTGGAGCAGAGGTCTGAAGCCCTTGACCCCTGGGGGCC
TGGGTAGTGTAGGATCTCGCTGGGCTGGTCTTCCATTCTTCCAACAGTTTCATGTTCACTACTGGACTTTAC
TATGGGATAAGGCCCCCTGGGGTCTCCATTTCTTCCAACAGTTTCATGTTCACTACTGGACTTTAC
GGGCTCAGTATCTCCCTTAGCCATGAGCTGGCTCAGGCATCCCTTCCCTTCCCTGGAGCTGCCCTGC
CTTTCTCAAGTATTTATTTATTTATTTGTCATGGTCTTGGGAACATGTGGCACAAGTAATGGGATGAGGA
GGAATTTGGGGTGGGGTCTTCTACCTAGGACTCTTCCCTGGAGTCATGGGCTGCCTGGGACCCAGGAC
CCATGAGGGGGCTGAGAGGTTTCTACACTCGAGGAGCAGGGTCCAGAGAGGCAGGCTGGGGAGGCAAG
GGACCCATCCTAGGCCGCTTTCTTGGCAGCCAAGCAGCTTAGCTGGGGCTGTGCAGCCAGGGGCTTA
CCCAGGCCAGTGGAGTGCCACAGCCCTGGGGAGCCAGAGGCTTTGGTATCGTATCGCCTCTGTGTC
CTTTAAGAGAGGAGAGTTCAGTACCCCGTCTTTCTTTACTGAGAGGAACTAAAAGGATCTCTGT
GTCTATGGAGAATTGTCAATAAAAAGGCCTCAAGCTTC
ACGCGTAAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_144588.7](#)

Summary: This gene encodes a protein with several transmembrane domains, a Rab11-binding domain and a lipid-binding FYVE finger domain. The encoded protein appears to promote neurite formation. A mutation in this gene has been reported to be associated with hereditary spastic paraplegia, however the pathogenicity of the mutation, which may simply represent a polymorphism, is unclear. [provided by RefSeq, Mar 2010]

Locus ID: 118813

MW: 57.9