

Product datasheet for **SC215743**

NBR1 (NM_031862) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	NBR1 (NM_031862) Human 3' UTR Clone
Symbol:	NBR1
Synonyms:	1A1-3B; IAI3B; M17S2; MIG19
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_031862
Insert Size:	1644 bp



[View online »](#)

Insert Sequence: >SC215743 3'UTR clone of NM_031862
The sequence shown below is from the reference sequence of NM_031862. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AACCAACGACTGGTACAGCCAACGCTATAGGAGTGACCTTGATTAAATAACTGCCTGCTGCTCAGAG
ATGATCTTTATTCTGTCAATTGGGGTATGGGATAGAAGCCCTTGCTATTTTTAATCTGATGAATCTGTA
TAGAGCCCATCGTTGAGTTACCAAGACAATACCTGCTACAGTATTTGGGGAGCAAATAAGACCAGA
ACTTAAATTTTCACTTTAGACATTGGATGAATAGTATGAAGACAGTTTTTCAGTTGATTTGGATAAAAC
TATTTTAGTGATTGACAAGTGAACCTCACTTCATATAGAACCATTTTCTTCTGCTTTTATTGAA
ACTGAGTATTTTCTTGGCTAATGTGGATTTTTATGGGGATATCTGTTAATTTTCAGTTTGTAAAG
ACATTAACCTCGGAAGTTGTTTTAAGAATTATCTCATAATTCTTATTCTCATAATTTCTGTAATCCA
CCTCAAGCTTCATAGTTATTTGGCATTGAAATAACACCCAGAGCATGATAGAAATGTTGTTACTCTTCC
TCTCTCAAGGAGAAAGTAATTTTCTGCAATACTTAATAATTGGCACCGTTGCTTCTAAAGACTCCAT
GGTGCAATCAAGAGTATCCAACCTCAAGGAATCTCCGCAATTTCAATGAAGGAGGAAGAGTGTGCTGA
TAAACCTACCAGCACCTATTGAGCAATGTCTATTATAGTAATTTTGCATACATTTTTATTTAAGGGAAA
AAATATAGGTATTGTGAAATATTTTGTCAATCTTAGAAAAGGAAAAATCCCCTTATTTAAAGGGAA
AAGTAAATTTAACAGTTGCCTTTTTCTTAATGTCAGGGCAGATCTTATTTACAGTACAGTGGGGGAA
ATAGAAACATGTGAAAGGCAAAAGGCAGGCTCCTAAATTAATGTCAGTGAAGTTCAGGGTGGGCAATG
AGTGTGTGTGAGGTATAGGAAATGCTGATGACTTCTTAAATGCTTGAAGTCCGTTACAGGTATCTAGC
CCTAGAATGCCTAGAACAGGAAGAGGCAGCTGGTGTCTGCAAACTTGGACAGGGGCAAGTTGCTGA
AAAAGTTTTGGTTTAACCCGAAGATAAGTGGAAGAGAGCTTGCCATGAACCCAGGTTCTCACTCTGTT
TACAGAAGTGTGTGAGTACAGTTGGTGAAGGAAGAGGTAACAAAAAATGCTAAATTTTTATCCATGA
AAATGACTTCCAGAAAAGGAAGAATATGAACCCAGACCGAAGGGGAAAAGATAGTTAATAGTATTATC
TAACCTGGTTGGTATTTGTAATGAATGGTGATTTTAATTAGTCATTAGCCATAATGATGTTTATTACA
GTATACTCCTGAATGCTACTTAAATAAACCAGGATTCAAAGTCAAGCCAGCCAGGCCGTTTATTATT
TAAACGTTTTAATCGGGGCTTCCGGGTAGAAGGTGGAGCGGCAGGGTGAATTGGGTTGATGGGTGGG
ACCTGTCTTGACCATCGGAGTTTTATAATCGAGGGCCAGGAGGGCCGGGTTGCTCTCCTGGTTATGTA
TGTAATTTGTACATAACCACCTAAAGAATGGTGAAATAAATGTTCTTGAAATTCCTA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: NM_031862.4

Summary:	The protein encoded by this gene was originally identified as an ovarian tumor antigen monitored in ovarian cancer. The encoded protein contains a B-box/coiled-coil motif, which is present in many genes with transformation potential. It functions as a specific autophagy receptor for the selective autophagic degradation of peroxisomes by forming intracellular inclusions with ubiquitylated autophagic substrates. This gene is located on a region of chromosome 17q21.1 that is in close proximity to the BRCA1 tumor suppressor gene. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Apr 2014]
Locus ID:	4077
MW:	62