

Product datasheet for **SC215761**

ZNF169 (NM_194320) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF169 (NM_194320) Human 3' UTR Clone
Symbol:	ZNF169
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_194320
Insert Size:	489 bp
Insert Sequence:	>SC215761 3'UTR clone of NM_194320 The sequence shown below is from the reference sequence of NM_194320. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC CATCAGCTCCTACCCCAAGAGGTCTTCTGACCTTTCTTTCCCGTGAGTGTGAAGCTGGCAGAAATCA CTAGTAAATGCTTCAGTTTCCTGAAATGGAATGGAAAAAATGTTGCTTTCTTTCATCGATCATG AGATAGTTGATTTTGGTTTACTTTCTATATTCACAATTTGTCTTGGCTTGCTAGGGGTTCCATAACAA AGTACCCAGGCTGGGTGATTTTGACAACGGAAATATATTTTCATATTTATGGAGCTGGAAGTCCCAG ATGAACGTATTGGTAGGTTTGGCTGCTTAAGCCTCTCTGTTTGGCTTGGGATGGCTGCTTTCTGG CCACGTCCTCACGTGGGCTTTTTCTGTGTGTATGCTTTCTGGTGTGGCTTCTCCTCTTCTAGGGAC ACCAGTCATTTTGGATAGAGCCACAGCCACATGGCCTCGTTTAACCTTAATAACCTTCTAAAGGCCTT ATCTCT ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.


[View online »](#)

RefSeq: NM_194320.4

Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Locus ID: 169841

MW: 17.9