

Product datasheet for **SC215883**

ATG4B (NM_178326) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ATG4B (NM_178326) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	ATG4B
Synonyms:	APG4B; AUTL1
ACCN:	NM_178326
Insert Size:	1682 bp



[View online »](#)

Insert Sequence: >SC215883 3'UTR clone of NM_178326
The sequence shown below is from the reference sequence of NM_178326. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGCTGCCAAGTCCAGATTCTTCTGATGTAGAGCGACTGGAAAGATTCTTCGACTCAGAAGATGAAGACT
TTGAAATCCTGTCCCTTTGAAAATCCTGGGGTCCGGGGTGGCACCTGTGAGAGCCTGGGGCTCCTGGTG
CCGCTGCGTTTTCATCCATCCCGCCGCTCGCCTGCCGAGGGCTGCGCCCGTGTGCTCCCGCCAGAG
GGCCACCCGCTGTGCTCGTGGACTGAGGCTGCGCTGCCGGGAGGCCCTTACTGCTGGTGTGACTGCTG
CCAGCTCAGAGTGCCCGTCAGGGCTGTGCATCCGCACGCGGAGCCGTCTGTTAGGAGCTCCAGAGTG
TTCTCTCGACACTGCCAGCCCGTGTAGCACCTGGGCCTCAGTCCACTTGCTCCAGGCGCCGGTTC
TGTGGTTGGTTTGAATTAAGTCTGTTGAAGTTGTCAGACACAGACATGAATTTCTGGGCGCTCCC
TGAGTCAGAGTCTCAGAAGACCTGTGCAGGCTGGCGTGAGAGGAGCGGCAGCCACTGCGGCCACG
CCCAAGGACTGGGCTGTCTCGAGGGGGCGCGCCACCCTGTGTCTCTGCCCAGCCTGGCTTAC
CAAGGGCTACCTCAGTGGGAGATGAGGTTGGAGGAACGAAGGCGAGGTTCTCTCTGCTTTGGGGAGAA
AAGTATTCAGGAAGTGGGTGTGTGGAAACCTGAAGATGGCGTGCACAGGACACAGCGTGGGCGGCCTG
GGCAGAAGGGCGGCTGGCTGTCTGGAGCTGTGCTGGAGCCTGCCCTCAGAGTGTCCCTTTCCAGTGC
TGTGGCATTCTGTGGCAGCTTCCCAGGTGTGGTGACGGGGGGGGGGCGGGGCTCCACCTGTGACAGC
CAGGCTTGAGGGTGGACGGCGTGCCTCTCCAGGAGCCTTCCCATGTCTTGCCTTGTGAGAATTGC
CCTCCCATGCCGCTGAGGTGTTAGTGGTTTAGGGCCAAAAGGGGAAAACCACTTGAGTCTTGTGGTGT
GTGGTGGGCAGACACCACAGGGTGGCATCACCTGGTGGCATTTCAGAACCTCAGCCCCGATTCCAGCA
CCACCACCGCCTGACCCTGTGTAACTGCTGTCCCGGGTCCCAGAGTGCACCTGCCCCGCTGTCTG
CTGCTGTCTGGGAAAGTATCTTTGCCCACTAGGAAATGTAACAGGAGGGCTTGGGGAGCGTGGGC
ACTTTTCTCATGAGCAGCTACTGCGCGTGGCAGGACTCGTGTGCTGCTGCTGCTTGTGTAGGTGCG
GGGAGCCAGAGATCCCCGAGGACGCGCGCCGACAGTCCGCACTGACCGCCACCTGGTAGCAGAGGA
CACCCCCAGCCCCAAGCATTGAAGACATAGTGTATTTCTCGTATCCTTTCTCCCTTGGGTGTAGTT
GGGGTGGGGAAGCAGGAAGGCTGGTGCATCTCCATTCTTGGGCTCCACGTCCGAGTTCATGGTGGC
CCGCTGTGCTGGGAGCTGCAGTGGTAATGTGTGGGACACCTTGACCAAAGGGAGCTTTGTCTCGTGTG
TTTTGAAAAGGCTTAATGAAGAGAATGTTGTTTATTCTTAGTAGTATAGTTGCAATTCTAATGGCA
AATAATAAGTTTCAGTAGAAAACAAA
ACGCGTAAAGCGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_178326.3](#)

Summary: Autophagy is the process by which endogenous proteins and damaged organelles are destroyed intracellularly. Autophagy is postulated to be essential for cell homeostasis and cell remodeling during differentiation, metamorphosis, non-apoptotic cell death, and aging. Reduced levels of autophagy have been described in some malignant tumors, and a role for autophagy in controlling the unregulated cell growth linked to cancer has been proposed. This gene encodes a member of the autophagin protein family. The encoded protein is also designated as a member of the C-54 family of cysteine proteases. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]

Locus ID: 23192

MW: 60.3