

## Product datasheet for **SC215934**

### **C15orf29 (KATNBL1) (NM\_024713) Human 3' UTR Clone**

#### **Product data:**

|                           |                                                   |
|---------------------------|---------------------------------------------------|
| Product Type:             | 3' UTR Clones                                     |
| Product Name:             | C15orf29 (KATNBL1) (NM_024713) Human 3' UTR Clone |
| Symbol:                   | C15orf29                                          |
| Synonyms:                 | C15orf29                                          |
| Mammalian Cell Selection: | Neomycin                                          |
| Vector:                   | pMirTarget (PS100062)                             |
| ACCN:                     | NM_024713                                         |
| Insert Size:              | 1695 bp                                           |



[View online »](#)

**Insert Sequence:** >SC215934 3'UTR clone of NM\_024713  
The sequence shown below is from the reference sequence of NM\_024713. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GTAGATGCTTATTATTACAGTTACATTGAGAGATTTTCATCTACTAAAGAGCATTGGTTTTTCAAAAC
ATCCCTGAACTGTATAATTTACAAAAAAGTCTCGTCTGAGAACTGTGAAGTGAAGAAATCA
AACTATTTTTCTTTAAAAAGCCACGTAATGAAACCACTAATGAAATCCCAGCAATCTGCTTCACAT
TGAAGTGGAAAAATATCCAAAAGGAGCAGCTTCAATTTTCATTGAGGTGAAAGTGCACATGAAGATTGT
TCACCTTTGCTGCATTTGGGAGTTATATGGTTATTTGGTAACATTAAGAACTACTGGATTTAATGCAA
TCCTGCATAAAAAATAATTTTACTATGTGAAAAATAAGACAGGACTTACCACTAGGAACCAACG
ACCAATCATCATTAACCTTTTTAAGATTGTGTTTTATTAACCACTTAAATGTGTGCAGCT
ATTTTCTTATGTGAAAAGACTGAAAGTTAAACATGAAAAAATCAATATTAACATTTTTTGTTC
CACTGAGATACTGTGTATGTAATGCCTTAATTATTAATAAGCCAATGTGTTATGATACCAATATCTG
TTTTAAAAAACTAAACCAACCATGCTTCTGGCATGATAAAATCATGGAATTAATCAGGGGTTACAT
TCTTGTAGAGTGTCTTGAAACACTCTCTGCACCATTTTTAAACCTTGAGAATAGTTTTAGTATCTCTG
ATATTTTTGCCAGAATCATCATGTATGATGAATGTGTTATCCCTATCTAAGGAAAAAGTGAATAT
GTTTTTGTATGAATGTTAACTGGAAATGTCCATGGACTTGGCTAATTTATTTACTTTTTATTGTAC
ATAGATTTCTAATATTTTCATTCCTGTATCATTTAACTTCCTTCATTTGAGTAAATCACTAAATAT
TTCTATTTTTGCTTTTTTAAATCTGATTTTATGAATCTAATCTTTTTCACTACATATGTTTTA
AAGAGTTACATACAGTGATTTAGAATGGTTTACAGTTAATGCTGATCTGTATTTTAAATCCAACACT
TTGTGTCACCTCCTCTAATGGTTAGTATGATATGCTAGCAGACTGTATGAGGTCTTTTTTAAAT
ACCACTTTTAGTGTGAGTGAACCAAACTTCTGGAATGTCTTAACAGCTCTAATCTTACTTGTCTTGA
ATGATTGGGGTTAATACCACTGCTGGTGGTTCACACATCATCCCATCTTAATATGCCTGACAGGCAT
CTGAGCAAAGGTTTTAGTAATTGAATTTCTCTGCAGTAGTCTTCAAGCACTTGAATGTAACCTTTA
GCATTTATTCGTTAATGACTACTGATACGAATCTCAAGCAGATTTCTTGTCTTAAAGTTATGTTTC
ACTGAGTTCTGGTTTGTGTAGCTATATTTATATAGCTAGATATTCCTCACAGTGAACATGAATTGTA
ATAATTGGTTATTTCTTAAGTCTTTAGATTATAATAATTCAGATTATTGCACGCTGTGATTTGAGA
GGTGAGTTATTTAAGAGGCCAGTTTCAGGACATGGGAATTTGAATTGTAACCTGTTATCTCTGTGAA
ACTTTTAACATGATAAAATATAACCTTTCTTTGTGCTTA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

**Restriction Sites:** SgfI-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_024713.3](#)

**Summary:** Regulates microtubule-severing activity of KATNAL1 in a concentration-dependent manner in vitro.[UniProtKB/Swiss-Prot Function]

**Locus ID:** 79768

MW: 67.3