

Product datasheet for **SC215953**

TIGD6 (NM_030953) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TIGD6 (NM_030953) Human 3' UTR Clone
Symbol:	TIGD6
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_030953
Insert Size:	1695 bp



[View online »](#)

Insert Sequence: >SC215953 3'UTR clone of NM_030953
The sequence shown below is from the reference sequence of NM_030953. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAAATTACAGATTTCCTCCAAACAAAAATATGCAGGAATTTATTTTCAGAAAATGTAGTTTACAAGAATA
AAGATTCTTTAGATAGTTTGTGAGCCAAATTAAGTAAAGCAATGTTATTGTGACAACATTCCAGTAC
TCTGAAATAGCCAGGAACTTCTTTGAATGGAATTTGACTAATATGTGTGTTTTCTTTCTTTGTTTT
TGGCTGTCTCTGGTCCTTGATTCAAGATGTATTTTGATTCATCAAGGGTTTCCAACTTGCTGCAAA
TTAGGATCACTTGAGAATCCTTTAAAAATTCCAAAGCTCAGGCCATATCCAGGCCTATTAATCACAA
TCTTTGGTGACGGGTACAGGCATTGGTAGTTTGAAGCTCTCCAGGTGATTCCAATGTGCAGACAAAT
TTGAAAACGAACCAACACAGGAACATCAAGTACACTGTGGCCTGGGTCCAGCTTCTTTCCAGTAAG
TGTATCTCAGGGCTTCCAAATTTAGCTTACATCAGAATCACTTGAAGGCTTGTTAAACCCCAAGGCT
GCTAGGCCACCCCGAGTTTGATACAGTAGACCTCAGGTGGGACCCAAAAATTTGCATTCTAACAC
ATTCTCAGCTGATGCAGTCCAGGCTCCATGCTTTGAAAACCACTGGTCTAGCTTTAGATAGGATATTGA
GCCAATTTAAATAAGCAATATACTGGTCTAGCTGAGTTCTAGAATTCTCTCTTTAGCTGGCCATCTG
AATACTCCCCATCACTAATTGTTAAAAAAGAATCAACTGTTCTTACTCTAGAGCTCTTTTCTCTTT
CTGCTGATTTGCTGGAAGCACTACAAGACTTCTGTTTGTTCGTTGTTTGTGTTGTTGTTTTAA
AGATGGGTCTTGTATATTGCCTAGGCTGGAATGCAGTGTTATTACAGGCATGATTATAACACACT
ACTCTCTCAACTCCTGGCCTCAAGCCATCCTCCAAATAGATGGGACTACTGATGCACACTGCCATGC
TGGCCTTACGAAATGTTTAAATAGGCATTTCACTAATAGGGATCTGGAGTACAAGGAAATACAGTGCAT
TTAAGACATAGGCTGGGCATGGTGGCTCATGCCTGTAATCCAGCACATTGGGAAGATCACTTGAGGCA
AGGAGTTTGAGACCAGACTGGCCAAACAGCGAGACCCCATCTCTAAAAAATAAATAAGAC
ATGGATACAACAGGGAGGGAGGTTATAATAATAGAGCATCCTGTTAATCAGAATCATGGGGAGAAACA
TACAGCCAGCGTGCCTCTCACTTTCCAGAAAGGGTGAGAGGATCATTAGAGAGTTTATTAGAGAGTAT
TTATTCCTAATAAGTCAGATAAATGTCAGTTTCATTTTTAGAGTTTACACAAATACCTGCTAATGGTG
TGATGTCATTTTTCCCCTGCTTCAGGTTAATTTTTGGATATTTTCAGAAACCCATAGCAATTCAGTGA
TTTTCTTTTGTAGTAGCTCAGTCAAAACAGGTAACAGCAGATTAATCATAGAAATGTTGCTTCCAC
ATTAACAAAACAATCACTGAAAAACAGCAGATGTTTAAATAGATTATTTATAGCTTATTTTGATATT
TTTTCTTTAATAAACTAATATGATCATTTGAATTAATA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_030953.4](#)

Summary:	The protein encoded by this gene belongs to the tigger subfamily of the pogo superfamily of DNA-mediated transposons in humans. These proteins are related to DNA transposons found in fungi and nematodes, and more distantly to the Tc1 and mariner transposases. They are also very similar to the major mammalian centromere protein B. [provided by RefSeq, Oct 2009]
Locus ID:	81789
MW:	65.6