

Product datasheet for **SC215969**

HNT (NTM) (NM_016522) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	HNT (NTM) (NM_016522) Human 3' UTR Clone
Symbol:	HNT
Synonyms:	CEPU-1; HNT; IGLON2; NTRI
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_016522
Insert Size:	1706 bp



[View online »](#)

Insert Sequence: >SC215969 3'UTR clone of NM_016522
The sequence shown below is from the reference sequence of NM_016522. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTGGTCTTGACCTGCTTCTCAAATTTTGATGTGAGTGCCACTTCCCCACCCGGGAAAGGCTGCCGCCA
CCACCACCACCAACACAACAGCAATGGCAACACCGACAGCAACCAATCAGATATATACAAATGAAATTA
GAAGAAACACAGCCTCATGGGACAGAAATTTGAGGGAGGGGAACAAGAATACTTTGGGGGAAAAAAG
TTTTAAAAAGAAATTGAAAATTCCTTGCAGATATTTAGGTACAATGGAGTTTCTTTTCCAAACGG
GAAGAACACAGCACACCCGGCTTGGACCCACTGCAAGCTGCATCGTGCAACCTCTTTGGTGCCAGTGTG
GGCAAGGGCTCAGCCTCTCTGCCACAGAGTGCCCCACGTGGAACATTCTGGAGCTGGCCATCCCAAA
TTCAATCAGTCCATAGAGACGAACAGAATGAGACCTTCCGGCCCAAGCGTGGCGCTGCGGGCACTTTGG
TAGACTGTGCCACCACGGCGTGTGTTGTGAAACGTGAAATAAAAAAGAGCAAGAAAGAAAAAGGAAACAA
AATAAGACCGTCTGACAGCAACAACGGTCCCACAAACAAGTCACAAAAGATACCGTTAAACTTTTTTTT
TTTTTATCATTTTACTACATGAACATCATGGATAACAAGGGATTCTGATTCAATTATTAATTTTCATTAA
TTATATTTTCTGATATGAGTCTAGAACTTACTGCAAAAAACAAGACAAAATAAAAAATACAACAGAGA
AGGGTGAAGAGAAGTATGTTTGTAAACAGTAAAAATGAATAGTATACTTCTTAAGTTTACAAC
GGTGGACCACACACAGGCACTAATCACCTGGTGAGGATTTGGCATATCCACCAAAAAATGCATCCGAT
TTAACCAACATCTCCACCAGCGCTACGGACTCCTCCCAATTCTGCATCTCTTGACAGACAATACTATGC
TCTCTACACACTGTTTAGAAATGGAAAGGTGATCTGCACTGTATCTTGGGTTTGTGGCTATGCTTCCT
TTGATGACATATATTATACAGTATATATATACATATATTTTTTTTGTAGAGTTCTAGCCATTTTATTT
CTCCGCAGGGTCCTTTCTCAGACATTACTGCATGCTGTATATGGCGTTAGCTGTGTGTTGATCTTCTAA
AAGATGATAGAGTTTACTGGTAATTGTGTAATCAGCTCCTGCCTTTTTATTTTCTGGGTTATTTACAT
GTCAGAGACATTTATAAAAAAGTAAAGGATAAAAAAAAAAAAAACAATAACCGGGCGCAGCATCT
TTCCAGGTGTGGCTCTGCGCTAGGCCACCCAGGTGGCCCGCTTTCTCCATTCTGCCCTGTGGCGTTA
ACAGACAGCAAGCAGCTGAACAAGCAGTACCGTCAGTACCCACTTGCTTTAGCCCATATGGGAATCTC
TGATGCCTTTGTGACCACTGGAGAGTTTAATCCTCTTGGTTTATTTTATTTCCACCTTTGTGTGAGTG
TGTATGAAAGAGAAGAAATGCAATTTTAGGTAATCTTTTTTTTTTTTCCCTCCCTGAAAGTCTGG
GAACCTGAGAGTCTCGGGGAGGGTCTGGATGTGATGAGGGAAGGGGATTGGGGATCCGGGAGG
GTGGGGTTGTCTGACTTGACATTAAAAAGTGTCCATGTCCCTCTTCA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_016522.3](#)

Summary:	This gene encodes a member of the IgLON (LAMP, OBCAM, Ntm) family of immunoglobulin (Ig) domain-containing glycosylphosphatidylinositol (GPI)-anchored cell adhesion molecules. The encoded protein may promote neurite outgrowth and adhesion via a homophilic mechanism. This gene is closely linked to a related family member, opioid binding protein/cell adhesion molecule-like (OPCML), on chromosome 11. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2009]
Locus ID:	50863
MW:	65.4