

Product datasheet for **SC215974**

VSIG8 (C1orf204) (NM_001134233) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	VSIG8 (C1orf204) (NM_001134233) Human 3' UTR Clone
Symbol:	VSIG8
Synonyms:	FLJ20442; FLJ39187
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001134233
Insert Size:	1690 bp



[View online »](#)

Insert Sequence: >SC215974 3' UTR clone of NM_001134233
The sequence shown below is from the reference sequence of NM_001134233. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site
Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**CGATCGC**

TCTCACCTGACCACAAGCATCCTGGAGCAAT**TGA**GTTTCAGGCAGGAAGACTCATCTGAGCTGAGCTGCT
TTTCTTGAAGTTGTGCAAGAACGTCCTGTTGGCTGGGGCAAATTCAGTATTTCTTTCTGGCCACTGTG
CTGGAGCGTGAGGAAAGACTCTCCACCCTCACTCCCAGATCCTGCAGGAGGCGTGAGAGCGGATGGAGAT
GGTGAGGATGTCTATTAATATGGTGTTCCTGCAAAACGGTGCCCATCTCAACCATGGGTTGTCAACCAT
CCTCGGTGTCTACAAGTCAATTTACCTTTGAGGGCCCATGAACTCTGATTACCTGAAAGTCCTTAGGCC
AGAAAGAGAAACCCCTGTCCCCACCTCTCACTTTACCTATTTCTAGACCTTTTCTTACAGCAATGCCT
TGGCCCTCAGGAAGCATCAGATCTCAAAGTGCAGGGTCTGGGGCTGAAGCTGATGCTGGCATGGGTGACT
CAGAGTTCCTCTATTAGTTCAAATATGTAATATCTTTAGAGCCTAGGACTCAGGCTGGAGGCAATTCTCT
TGGTGCCAGTTGGTTAATACTTTATGAATAGTGGTGGGGGAGGGAATGCCTTGGACTATTAATGTCCCCA
TGGATCTGGATCCTGGAGAACTCCAGTGTCTGGATGGAATGTTCTGCTTGGCTCATCTTGCCAGAGA
AGAGCCAGGATATGAGAAGATGGGAGAGGGAAGGGGAGGTGTGAAACGTGCTGGTGCTTCTGTCTCTCG
TGGTGTGATGGCTGCCAGGCTGCTCCTTCGAGGCAGGACAGGATTCCCATGGGAGATATGAGGACCATG
GTGATTCCAGGCACCAAGTTATGTCCCAGAGCATCGTGTGCCTGTGCAGTGGGGTTACTGATAGTCCAGG
CATCAGTTCATATTGTCTTTATGGCAGATCCTGCACAAGGGGTTCTCTGTCTCTGGACCCACTCTGTCT
AGCATGGACATCCTTTTCTTTTCCCTGGTAAGGACATCAGAAAGAAAGTTACCCCTCCAGTGTGTTGGA
GTCTTGGAAATCCTTGGCCTGCCTGCAGGCCTCTGAATCCTCAGTGTGTTAGCCCAAGGTACCTTAGGTA
AAGGAACATTTCTCTGGTCTTCCAACACCCTGAAAGGGAGGCTCTTCAAGAATAGGAAATTTCTCCTA
GTATCTAACCTTGGCCTGCAACTATGCCATACATTTAATTGTAGAAGGTATGGAAGTTAGGACTGTGGT
GGACCCAGTTAAATCACTGGACAGAAGCCATGAGGCAGAGCCAGAACCAGCTTCTGGACTCCCTGCATGG
AGCTTTCTACTATCCCAGGTTGTTGGGTCCCTTCTCAGTGGCCAAAGATGCTATTTTATGAGCCTTGAGG
ACATGACACTAAGTGAAAAAGCCTGGCCAAAGGGCCACATAGTGTATAGTTCTATTTATGAAATGTC
CAGGAAAGACAAAAGCATAGAGGCAGAAAGCAGACTAGTGTGTGGGAGGAATGGGGATTGGGAGGTG
ATGGCTAAGGGTTATGGGGTTTGGGGGAGGTAAAGAAATGTTTTAAATGGTTGTGAAGGTGGATG
CCCAACTCTGCAATACACTAAAGGCATCGAATTGTATACTTTAAATAGGTGAATTGTATGCTATGTGA
ATTATATCTC

ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAATGACCG

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: NM_001134233.1

Locus ID: 284677