

Product datasheet for **SC216095**

alpha 1 Antitrypsin (SERPINA1) (NM_001127705) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	alpha 1 Antitrypsin (SERPINA1) (NM_001127705) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	SERPINA1
Synonyms:	A1A; A1AT; AAT; alpha1AT; nNIF; PI; PI1; PRO2275
ACCN:	NM_001127705
Insert Size:	1717 bp



[View online »](#)

Insert Sequence: >SC216095 3' UTR clone of NM_001127705
 The sequence shown below is from the reference sequence of NM_001127705. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site
Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**GCGATCGC**

CCTCTTCATGGGAAAAGTGGTGAATCCCACCCAAAAA**TAA**CTGCCTCTCGCTCCTCAACCCCTCCCCTCC
 ATCCCTGGCCCCCTCCCTGGATGACATTAAGAAGGGTTGAGCTGGTCCCTGCCTGCATGTGACTGTAA
 TCCCTCCCATGTTTTCTCTGAGTCTCCCTTTGCCTGCTGAGGCTGTATGTGGGCTCCAGGTAACAGTGCT
 GTCTTCGGGCCCCCTGAACTGTGTTTATGGAGCATCTGGCTGGGTAGGCACATGCTGGGCTTGAATCCAG
 GGGGACTGAATCCTCAGCTTACGGACCTGGGCCATCTGTTTCTGGAGGGCTCCAGTCTCCTTGTCT
 GTCTTGGAGTCCCAAGAAGGAATCACAGGGGAGGAACAGATACCAGCCATGACCCAGGCTCCACCAA
 GCATCTTCATGTCCCCCTGCTCATCCCCACTCCCCCACCAGAGTTGCTCATCCTGCCAGGGCTGGC
 TGTGCCACCCCAAGGCTGCCCTCCTGGGGGCCCAGAACTGCCTGATCGTGCCGTGGCCAGTTTTGTG
 GCATCTGCAGCAACACAAGAGAGAGGACAATGTCTCCTCTTGACCCGCTGTACCTAACCAGACTCGGG
 CCCTGCACCTCTCAGGCACTTCTGGAAAATGACTGAGGCAGATTCTTCTGAAGCCATTCTCCATGGGG
 CAACAAGGACACCTATTCTGTCTTGTCTTCCATCGCTGCCCGAAAGCCTCACATATCTCCGTTTAG
 AATCAGGTCCCTTCTCCCCAGATGAAGAGGAGGGTCTCTGCTTTGTTTTCTATCTCCTCCTCAGACTT
 GACCAGGCCAGCAGGCCCCAGAAGACCATTACCCTATATCCCTTCTCCTCCCTAGTCACATGGCCATAG
 GCCTGCTGATGGCTCAGGAAGGCCATTGCAAGGACTCCTCAGCTATGGGAGAGGAAGCACATCACCCATT
 GACCCCGCAACCCCTCCCTTCTCCTCTGAGTCCCGACTGGGGCCACATGCAGCCTGACTTCTTTGTG
 CCTGTTGCTGTCCCTGCAGTCTTCAGAGGGCCACCGCAGCTCCAGTGCCACGGCAGGAGGCTGTTCTGA
 ATAGCCCTGTGGTAAGGGCCAGGAGATCCTTCCATCCTCCAAGGCCCTGCTAAAGGACACAGCAGCCA
 GGAAGTCCCCTGGGCCCTAGCTGAAGGACAGCCTGCTCCCTCCGTCTTACCAGGAATGGCCTTGTCT
 ATGGAAGGCACTGCCCATCCCAAATAATCTAGGAATCACTGTCTAACCCTCACTGTCATGAATGTGT
 ACTTAAAGGATGAGGTTGAGTCATACCAAATAGTGATTCGATAGTTCAAATGGTGAAATTAGCAATTC
 TACATGATTCAGTCTAATCAATGGATACCGACTGTTCCACACAAGTCTCCTGTTCTCTAAGCTTACT
 CACTGACAGCCTTCACTCTCCACAAATACATTAAGATATGGCCATCACCAAGCCCCCTAGGATGACAC
 CAGACCTGAGAGTCTGAAGACCTGGATCCAAGTTCTGACTTTTCCCCCTGACAGCTGTGTGACCTTCGTG
 AAGTCGCCAAACCTCTCTGAGCCCCAGTCATTGCTAGTAAGACCTGCCTTTGAGTTGGTATGATGTTCAA
 GTTAGATAACAAAATGTTTATACCCATTAGAACAGAG

ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCG

Restriction Sites: Sgfl-Mlul

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001127705.1](#)

Summary:

The protein encoded by this gene is a serine protease inhibitor belonging to the serpin superfamily whose targets include elastase, plasmin, thrombin, trypsin, chymotrypsin, and plasminogen activator. This protein is produced in the liver, the bone marrow, by lymphocytic and monocytic cells in lymphoid tissue, and by the Paneth cells of the gut. Defects in this gene are associated with chronic obstructive pulmonary disease, emphysema, and chronic liver disease. Several transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Aug 2020]

Locus ID:

5265