

Product datasheet for **SC216135**

YPEL5 (NM_001127399) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	YPEL5 (NM_001127399) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	YPEL5
Synonyms:	CGI-127
ACCN:	NM_001127399
Insert Size:	1720 bp



[View online »](#)

Insert Sequence:	<p>>SC216135 3'UTR clone of NM_001127399 The sequence shown below is from the reference sequence of NM_001127399. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site</p> <pre>GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GAGGAGCATGTACCATCTGATAACTCTGAAGATACAGAGAGAAATCCATCTTTTCCAGGTCTCCTTC ACTGAAAACAAAAATCTACTTACATACACTGTACCTTAGCATCAGAGTCGGATTAATGAACTGCGGAA CAAGAGGTTGTGAGAATCTAAGATGGAACCTTTCTTTCTTTCTTTTCTTTTAAATTTGTATTTT CCATCCAACAGCAGTGTGTAGAGAGAATATTATGCAGATGCCGTTAATTTTTACCCTATGTTTACATC TTGAGGCAGCAGAGTCTGTCTGCAGCTATGTGGTGAGCTATGTAAGGAAAAAATCTGGGCTGTTAGAG TGAAAAAGTGTGTTTTATGTCAATTGTGAAAGGAAAAATGTTAGGAGTATGGTTTTAACTGGGCTTC ATTTTAACTTTTTTTTTAAACCCAGTATTTCACTTGATTTGCTAGCTTCAGAGAAGAGATCCGAAT CTGTGCCAGCGCTAAAGGCTCAGTGTTAGCATGGCTTGTGCTGGCCGGTGTCCATATTCTTGTGGA GATGAACCGTAGCACCAGAGCCATTCTCCTTGTCAAGTCTGGCCAAAGATGTACCATTCTAGTT ATTTGTCACCACATAATTGGTGTGATTGGAACTTTTTCTGAGATGGGACAGAAGTCTGGGTTGTCT TTTTCCATGTAACCTAAGCATAGTAATATAAAATAAGTAATAGTTGGATGCTTTTGGTCTGTGTTGCT TTTAAAAACACCTTATAAAAGAGGAGAGATTTTGATAAGCAATTTTCATAGTAGTAAAGTTTTTTTCA TCTCTTAACTAAATTGACCATGCATATAATATCTTTGTTTAAATGAAAGCATACTGTTGAAACCCGC AGTGTTCATTTAGAAAACAGTTGAACAGAATGTCAATGTGCATTCATGCAAAAAACATTTAATCTGC ATCTGTTTTAGAAAAGGGGAAATGAAGCAACTGTCTAAAAATACTGCTTTACAAGCATTTCAGCCT TTCCCCCTCAGTTTTGCATTGATTTTTGACAAGTCTGTAGAGCCTAATAGTTTCCATCAAAGGCCTAG ATCTCTTATTTAGCATTTTTTTTCAGCTCTTCTCTCAGAAGTTCAGCTGTTGAAACGAAAACTGTACTTT GTACCCTCACATACAAAGGGATCAAAATTTGACCTGGTGTATTTTAGCCCAAATTTATGACATTACAC AATATTAATAATGTAATGTTTCTTTACCCAACTACTTCTAGATATTCTAGTATTTGCTTCTGGTGGAA TTAAATGACGGTAAAATTGGCTAATTTTGAATGAATGAATGGATGGATGTTTTGCATGCTCAATTTT TAGGTCCTTTGTCTAGAAAGGAAATTTGCCTCAGTTGAATTAGTGAAATATTTCTGTCGTTGATTTAA AAGTGACTTCTGAGTACAGTTAAGTTCCTCCTATTTGCCACTGGGCTGTTGGTTAGAAGCATAGGTAAC TGATTAAGTAGGTATGATACTGCATTTGAAATAAGTGGACACAACTATCCTTTTCCACCATGGACTC AATCTGAGAACAACAGCATTTCATTTCCATTTCATTTCCATACTGGCTTTTGATTATATGCAGATTCCTAG TAGCATGCCTTACCTACAGCACTATGTGATTTGCTGTACAATAAAGTATATTTTGTCTTGCA ACGCGTAAAGCGGCCGCGCATCTAGATTCGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG</pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001127399.2
Summary:	Component of the CTLH E3 ubiquitin-protein ligase complex that selectively accepts ubiquitin from UBE2H and mediates ubiquitination and subsequent proteasomal degradation of the transcription factor HBP1 (PubMed:29911972). Required for normal cell proliferation (By similarity).[UniProtKB/Swiss-Prot Function]

Locus ID: 51646

MW: 66.5