

Product datasheet for **SC216151**

UNC5B (NM_170744) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	UNC5B (NM_170744) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	UNC5B
Synonyms:	p53RDL1; UNC5H2
ACCN:	NM_170744
Insert Size:	2000 bp



[View online »](#)

Summary: This gene encodes a member of the netrin family of receptors. This particular protein mediates the repulsive effect of netrin-1 and is a vascular netrin receptor. This encoded protein is also in a group of proteins called dependence receptors (DpRs) which are involved in pro- and anti-apoptotic processes. Many DpRs are involved in embryogenesis and in cancer progression. Two alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Oct 2011]

Locus ID: 219699

MW: 73