

Product datasheet for **SC216380**

C16orf55 (SPATA33) (NM_153025) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C16orf55 (SPATA33) (NM_153025) Human 3' UTR Clone
Symbol:	C16orf55
Synonyms:	C16orf55
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_153025
Insert Size:	991 bp
Insert Sequence:	>SC216380 3'UTR clone of NM_153025 The sequence shown below is from the reference sequence of NM_153025. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GACGCCTATAATTACATCTCAAAGAAACAAGCACCTTTGCATGGCACTCGCTACTCCCTGCGATA
GGCGTCTCGGGAACAGCCGCTCACCCTGTGAGAAGCCGAGGCCCTTCTCCAGTGCTCTCGGGAGGT
TGCACCAAGGCCACCCACCTGTGAGAACTGCAGCCCCTTCTCCAGTGCTTTCGGGAGGGTGCAC
CAGGCCCCCCCCACCTTGTGAGAAGCCGTGGCCCCCTTCTCCAGTGCTCTCAGGAGGGTGCACCAAGC
CTGCCCCCGCCGTGAGAACTGCAGTCCCCTTCTCCAGTGCTCTCGGGAGGGTGCACCAAGGCCACCC
CACCTGTGAGAACCTGCAGCCCCCTTCTCCAATGCTCTCGGGAGGGTACACCAAGGGTGTCCCCACCC
TGTGAGAAGCCACCGCCCCCTTCTCCAGTGCTCTCGGGAGGGTGCACCAAGGGTGTCCCCACGCTGTAA
GAAGCCGCGCCCCCTTCTCCAGTGCTCTCCAGCCGCCCCACCTGTGAGAAGCCGTGGCCCCCTTCTCCAGTG
CTTTCGGGAGGGTACACCAAGGCCCGCCCCACCTGTGAGAAGCCGTGGCCCCCTTCTCCAGTGCTCTC
AGGGAGGGTGCACCAAGCCTGCCCCGCTTGAAGAACTGCAGCCCCCTTCTCCAGTGCTCTCGGGAGGGTGCAC
GGTGCACCAAGCCTGCCCCACCTGTGAGAACTGCAGCCCCCTTCTCCAGTGCTCTCGGGAGGGTGCAC
CAGGCCCCCCCCACCTGTGAGAAGCCGCGAGCCGTTCTCCATTGCCATTGGTCCCTGCCCTCCCCAGA
GCTGCTAGCAGCTCGGTGAGGCAAGAGTGACTGGCTCAGAGAGCGAGCTTTGTGGTATTGATAAATGA
ATAAATTAGTATGTTTATGTTGGAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_153025.3</u>
Locus ID:	124045
MW:	34.3