

Product datasheet for **SC216554**

PCDHGC3 (NM_002588) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PCDHGC3 (NM_002588) Human 3' UTR Clone
Symbol:	PCDHGC3
Synonyms:	PC43; PCDH-GAMMA-C3; PCDH2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_002588
Insert Size:	1832 bp



[View online »](#)

Insert Sequence: >SC216554 3'UTR clone of NM_002588
The sequence shown below is from the reference sequence of NM_002588. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAGAAAGTCGGGCAAGAAGGAGAAGAAGTAAACATGGAGGCCAGGCCAAGAGCCACAGGGCGGCCTCTCCC
CAACCAGCCCAGCTTCTCCTTACCTGCACCCAGGCCTCAGAGTTTCAGGGCTAACCCCCAGAATACTGG
TAGGGGCCAAGGCCATGCTCCCTTGGGAAACAGAAACAAGTGCCAGTCAGCACCTACCCCTTCCCCC
CCAGGGGGTTGAATATGCAAAAGCAGTTCCGCTGGGAACCCCATCCAATCAACTGCTGTACCCATGGG
GGTAGTGGGGTTACTGTAGACACCAAGAACCATTGCCACACCCGTTTGTGTACAGCTGAACCTCTCC
ATCTTCCAATCAATCAGGCCCATCCATCCCATGCCTCCCTCTCCCAACCCACTCCAACAGTTCCTC
TTTCCCGAGTAAGGTGGTTGGGGTGTGAAGTACCAAGTAACCTACAAGCCTCCTAGTTCTGAAAAGTT
GGAAGGGCATCATGACCTCTTGGCTCTCCTTTGATTCTCAATCTTCCCCAAAGCATGGTTTGGTGCC
AGCCCCCTCACCTCCTTCAGAGCCCAAGATCAATGCTCAAGTTTGGAGGACATGATCACCATCCCCA
TGGTACTGATGCTTGTGGATTTAGGGAGGGCATTGTTGCTACCAAGCCTCTTCCCAACGCCTGGGGAC
CAGTCTTCTGTTTTGTTTTTATTGTTTGACGTTTCCACTGCATGCCTTGACTTCCCCACCTCCTCCT
CAAACAAGAGACTCCACTGCATGTTCCAAGACAGTATGGGGTGGTAAGATAAGGAAGGGAAGTGTGTGG
ATGTGGATGGTGGGGGCATGGACAAAGCTTGACACATCAAGTTATCAAGGCCTTGAGGAGGCTCTGTA
TGTCCTCAGGGGACTGACAACATCCTCCAGATTCCAGCCATAAACCAATAACTAGGCTGGACCCTTCCC
ACTACATAATAGGGCTCAGCCCAGGCAGCCAGCTTTGGGCTGAGCTAACAGGACCAATGGATTAACTG
GCATTTCAAGTCCAAGGAAGCTCGAAGCAGGTTTAGGACCAGGTCCCTTGAGAGGTCAGAGGGGCCTT
GTGGGTGCTGGGTACTCCAGAGGTGCCACTGGTGGAAAGGTCAGCGGAGCCCCAGCAGGAAGGGTGGGC
CAGCCAGGCCATTCTTAGTCCCTGGGTTGGGGAGGCAGGGAGCTAGGGCAGGGACCAATGAACAGAAA
GTCTCAGCCCAGGATGGGGCTTCTTCAACAGGGCCCTGCCCTCCTGAAGCCTCAGTCTTCACCTTGC
CAGGTGCCGTTTCTTCCGTGAAGGCCACTGCCAGGTCCCCAGTGCGCCCCCTAGTGGCCATAGCCT
GGTTAAAGTTCCCCAGTGCCTCCTTGTGCATAGACCTTCTTCCCAACCCCTTGTGCCCTGGGTCCC
CGGCCATCCAGCGGGGCTGCCAGAGAACCCAGACCTGCCCTTACAGTAGTGTAGCGCCCCCTCCTCT
TTCGGCTGGTGTAGAATAGCCAGTAGTGTAGTGGGTGTGCTTTTACGTGATGGCGGGTGGGCAGCGGG
CGGCGGGCTCCGCGCAGCCGTCTGTCCTTGATCTGCCCGCGGCGCCCGTGTGTGTTTTGTGCTGTGT
CCACGCGCTAAGGCGACCCCTCCCCGTACTGACTTCTCTATAAGCGCTTCTCTTCGCATAGTACG
TAGCTCCCAACCCACCTCTTCTGTGTCTCACGCAAGTTTTATACTCTAATATTATATGGCTTTTTT
TCTTCGACAAAAAATAATAAACGTTTCTTCTGAAAA
ACGCGTAAGCGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCGCCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_002588.4](#)

Summary:

This gene is a member of the protocadherin gamma gene cluster, one of three related clusters tandemly linked on chromosome five. These gene clusters have an immunoglobulin-like organization, suggesting that a novel mechanism may be involved in their regulation and expression. The gamma gene cluster includes 22 genes divided into 3 subfamilies. Subfamily A contains 12 genes, subfamily B contains 7 genes and 2 pseudogenes, and the more distantly related subfamily C contains 3 genes. The tandem array of 22 large, variable region exons are followed by a constant region, containing 3 exons shared by all genes in the cluster. Each variable region exon encodes the extracellular region, which includes 6 cadherin ectodomains and a transmembrane region. The constant region exons encode the common cytoplasmic region. These neural cadherin-like cell adhesion proteins most likely play a critical role in the establishment and function of specific cell-cell connections in the brain. Alternative splicing has been described for the gamma cluster genes. [provided by RefSeq, Jul 2008]

Locus ID:

5098

MW:

65.4