

Product datasheet for **SC216560**

RPS24 (NM_001142285) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RPS24 (NM_001142285) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	RPS24
Synonyms:	DBA3; eS24; S24
ACCN:	NM_001142285
Insert Size:	1826 bp



[View online »](#)

Insert Sequence: >SC216560 3'UTR clone of NM_001142285
 The sequence shown below is from the reference sequence of NM_001142285. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CCTTCCTGCCTCTCTGGTCTGGGCCACTAATGTCACTGCCATGGCCGCTTGCTGCATTTCTGAGGATG
CTTCATCTCTCCACCTTCTTCTCCACTCAGCAGCCAGCAGGGCACTGTGGAAATCGGAGTCACATGAGC
TGGCACCTCTGTTCAGAACCCTCCAGGGCTCCACATCTCTCTCACCCAAATGCCAAAGACCTCCCCACG
CCCCACAATCCCCACGACCTGGCCACTGGCCTCCCACCACCTTCCAGCTCCAGCGGCTCTACCACA
TTTAAGGCTTTCCTTCTAGTTTAAATTTTTCTCGTCAGCAGTTGATTTTATTATTTCTGTATTATT
GGTATTTTCCACTAGAAATGAAGCTGCGTGAAGTTAGAGATTTTTTTTTTTGGTCTGTGTTCTTAATT
AGCTCATTGCTATACCCCTGGCGCCAGAACAAATGCCTTGGACACAGTACGCAGTAGACTAAATAAATA
CTTGTTGAATGACTGACTGACGGAATGACGGCTGTGTGGGGAGTGGATTGGTCTGTGAGGCAGAGGCTG
CGGTGGAAACTCAGGCAGGAGGTGATGGTGGTTCTTGGGGCTGCGGAATGCCAAGTTTAGAAGCTCTTC
CTCTGCTGTGGCACATGAACCGGTCACTCGAGAAGGCTTTTAGATTTACTTTGCCTAATCCCCTCTTAG
TGCATGTGGGGAAACTGAGGTACACAAAAGGAATCCCCACCAAGTTAGGGGCAGAACCTAGCCCCCTT
GTCTCCAGATGGATATCTTCTTTTTTTTTGAGACGGAGTCTTGCTCTGTTGCCAGGCTGGAGTGCA
GTGGTACCATCTTGGCTCACTGCAACCTCTGCTTCCCAGGTTCAAGCGATTCTCCTGCCTCAGCCTCCT
GAGTGTCTGCGATTACAGGTGCACACAACCACGCTGGCTAATTTTTGTATTTTAGTAGAGACGGGGT
TTCACCGTGTGGTTCAGGGTGACCTCAAACCTCTGACCTCATGATCCACCAGCTCAGCCTCCCAACGT
GCTGGGATTACAGGCATGAGCCACCGTGGCTGGACATCTTGTTATTAAGCTTCTTCTCTCTTTG
TAGGGGAGGGGGAGATGCCTCTGGTGGAGAAGACCAGTGTGGCAGTGAAGTGTCTGTTAGTGAACCTG
GTGGCTGGTTGAGGGTCTGCTGTTGGTGAAGTACTGAGGACACATACAAAGTGGTTTTCTCAGTGGTCACTTG
GTGTTGGTGAATAAGGGTCAGAAGATGGCTCCTGTCTAGGGCACTGCCAGTCGGTTTGAAGCTGAAA
TGCCTGCTTAGCAGTTTGAAGAACACAGACCTTGGAGGATCTTCTGGTTGCCTTCAAGAATTCATT
CTATCCCCTTCTGCTCCCCAAATTTGCTTTTCTTGGGGTGGTCTTGGTTGGCTAAGCCAAGAAAGT
ATGGCATCTACTCTCCATAGCAATAGCTCAGGAATAGGCAGTGACCCAGACCTGAACCAATCAGTGC
ATGGAATTACCCCTGGCCAAAGTGGTTGATTGAGGCTGGGTGCAAGCAGAGTTGTGAGAAGGCTCCCAT
TTGGTGGTTGGAGAGATCGCACTTGTCCAGAGGTCATAATGTGCAGATCTGAGGCTTGGAACTGCTGC
AGACATTTTGTACCACAAGTGAAGCCACCTGACGACACAGTTGACAATTTGGAGCAGGGCAGAGCTG
AGAGAACAGCAGGAAACAGCCAGAGTCTTGTCAAGCCTCCCTGAAGTATCTATACCCCTGGACTCA
GTTATGGGGGCTAATAAATGTTATATACTGTT
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001142285.2](#)

Summary:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S24E family of ribosomal proteins. It is located in the cytoplasm. Multiple transcript variants encoding different isoforms have been found for this gene. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Mutations in this gene result in Diamond-Blackfan anemia. [provided by RefSeq, Nov 2008]

Locus ID:

6229

MW:

68.2