

Product datasheet for **SC216592**

BHLHB9 (NM_001142524) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	BHLHB9 (NM_001142524) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	BHLHB9
Synonyms:	GASP3; p60TRP
ACCN:	NM_001142524
Insert Size:	2000 bp



[View online »](#)

Insert Sequence: >SC216592 3'UTR clone of NM_001142524
 The sequence shown below is from the reference sequence of NM_001142524. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAAGTTAAAGAGATTATTGAAACAATGTAAATGAGCCAGAGATAGAACATTTTAAGCCATCTTCAAAC
TCTAGCAGGCTGTACATTACAGTGTACACATTATACACTGCATCTTTAACACAGAGTCACCTGTGACAG
GCTCTAGTGTGAGCTAGACTATTTTGGGGTATCAAATGAATATTATACCTTGGGCTGAAAATGTTTG
ATTTTTATCTGTCTAGATTGGCATATTTTAAACATTTTACTTAAAGATAGCAAACAGTTCGTTTTAAG
TAAGCTAACTTGTTCAATAGTATCTGTGGCTTAAATGGCAAAAAGAAAATATCCTTGAGTTTGTAAAT
CTAGTTACAGAAGTAAGGCATACACACACAAAGATAACAGTACCTAGAGAGAGAGTGTGTGTGAGTG
TGGTGTCTCTGTGTGCACGTGCACGCTCATGGCCAAATGTGCGCACTCTACATAAAGGAGGCAGGA
GTTCTATAGGCTATTTAATGTAAGAGAACTATTTTTCTCCTGTCCAGCTGTATCAGATACTCGTTC
CGCAACACAGAAATGACTCAGAATCTCAGACAAAATGTATTATTTGTTCAATTTTAAATTTGCTACTAC
ATTCATAACTCTTAAATTGTTAGGCTGTTTCATTTACATCAAAGTTATCTCACAAAAGAGAAGGCAGGA
AACGTTTTGTGAGTGCCTATTCTATGTCAAACACTGTGTTGGCACCATATTTTACAAGTTTTTTCTC
TTCTCACAGTGATCTTGTGAGTTAGTACTTATATTTTTATTAGAATCATTATTTCTGGGTACCCTCCA
ATGAGAATTAGAGAGGTTAAATACCTTTTCTAGATTCCACAGCAGGAAGGTGGCATAGCTGTTTTGT
CTGACACCAGAACCCTCACCACACTGCTTTACAGTCTTCTGAAGGACATTTTGGAGTGGGGGGGC
CTTCAAAGCTCAGAGACTGGTTGAATGGTTTAAATTTTGAATGGATCATGTCCATGCCAGGTGTTACA
ATTTCTTAACTTCTCAAATTCGTGTGCCATTAGACATTTGGCTACATCTGGCTGGAGTCCAGGAGAA
AATTTGAGGTAATAGATGGATTTTATCTGGCAGTGTGAAAATAGTAGGAGCCTAAAACCTTTGTTC
ATGAGTAAGATGTAGATTTGGAGTCATCAGCAGGGAACTTGTAGATGCCTTTGTATCTTTTTCTA
TCCATTTGTATTTTTCTTTCCCAAAGAGCTGAAATCATTTTATACCCTGAATTTAAAAGGTAATTGCT
TTTTCTCCTTGCAGAAATCATACCTTAAATTTTTTTTTTTTTTTTTTTTACCATTATGATTTATCTC
CTATGTGCATATGAAAGGCCAAGGGCTTTATCTTTACTAGTAGGCAAGCCTGGGCATTTCTGCATT
TGTGCCAGCTAACACTACAAAATAGACTTAATTAGAGAAATTTAGACCCAGAATTGCAGCTTCAACAAC
ACAAAACAGGAAAGCTAGAAGAGTTGCTGTGCAACCTACAAGATAGTGATTCTTGAAGTTGGTATG
GACCAGAAAGCATTCCCAGGACATTTAGGCTATGGCTGAGTTCTCTAAGACTCAGGTGGGAATTATAA
TCCTGGTCTAAGCTGATATGAAAGTTTCTATGTGTGAAAATATGGAGTTATATGCTTTGAGATTTCT
GCCAGTTAACTAAAACAAGATCATACTGTAAATGTTTGAACCTGCTTCCCCCACTGTATACTGTAA
ACATCTTTGCATGTCAATAAATATGCCTCTACAACATATTTTGAATCACTTAATATTTTCCCATGTT
AACAGTTGGTATATAGTGTGTAACCTTTATTGCATATTTCTATAACATACAATATACACATTCTTTT
GCAAAAATGGGAACACAGTATTGAGACTGCCTTTTCTCTTAACTTATCTTTTCAATTTAGTATT
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001142524.2](#)

Summary: This gene is a member of a gene family which encodes proteins with a basic helix-loop-helix domain. Other members of this gene family encode proteins which function as transcription factors, either enhancing or inhibiting transcription depending on the activity of other DNA binding proteins. The coding region of this gene is located entirely within the terminal exon. The encoded protein may be involved in the survival of neurons (PMID: 15034937). Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Sep 2011]

Locus ID: 80823

MW: 77.4