

Product datasheet for **SC216717**

PAI1 (SERPINE1) (NM_000602) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PAI1 (SERPINE1) (NM_000602) Human 3' UTR Clone
Symbol:	PAI1
Synonyms:	PAI; PAI-1; PAI1; PLANH1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_000602
Insert Size:	1835 bp



[View online »](#)

Insert Sequence: >SC216717 3'UTR clone of NM_000602
 The sequence shown below is from the reference sequence of NM_000602. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTTTTCATGGCCAAGTATGGAACCCTGACCCTGGGAAAGACGCCTTCATCTGGGACAAAACCTGGAG
ATGCATCGGAAAAGAAGAACTCCGAAGAAAAGAAATTTAGTGTTAATGACTCTTCTGAAGGAAGAGA
AGACATTTGCCTTTTGTAAAAGATGGTAAACCAGATCTGTCTCCAAGACCTTGGCCTCTCCTTGGAGG
ACCTTTAGGTCAAACCTAGTCTCCACCTGAGACCCTGGGAGAGAAGTTTGAAGCACAACTCCCTTA
AGGTCTCAAACCAGACGGTGACGCCTGCGGGACCATCTGGGACCTGCTTCCACCCGCTCTCTGCC
CACTCGGGTCTGCAGACCTGGTCCCCTGAGGCCCTTGCAGGATGGAACACGGGGCTTACAGGAGC
TTTTGTGTGCCTGGTAGAACTATTTCTGTTCCAGTCACATTGCCATCACTCTTGTACTGCCTGCCACC
GCGGAGGAGGCTGGTGACAGGCCAAAGGCCAGTGAAGAAACACCCTTTCATCTCAGAGTCCACTGTGG
CACTGGCCACCCCTCCCCAGTACAGGGGTGCTGCAGGTGGCAGAGTGAATGTCCCCCATCATGTGGCCC
AACTCTCCTGGCCTGGCCATCTCCTCCCCAGAAACAGTGTGCATGGGTTATTTGGAGTGTAGGTGAC
TTGTTTACTCATTGAAGCAGATTTCTGCTTCTTTTATTTTATAGGAATAGAGGAAGAAATGTCAGAT
GCGTGCCAGCTCTTACCCCCAATCTCTTGGTGGGAGGGGTGACCTAAATATTTATCATATCCTT
GCCCTTGAGTGTCTTGTAGAGAGAAAGAGAACTACTAAGGAAAATAATATTTAACTCGCTCCTAG
TGTTTCTTTGGTCTGTGTACCGTATCTCAGGAAGTCCAGCCACTTGACTGGCACACACCCCTCCGG
ACATCCAGCGTGACGGAGCCCACTGCCACCTTGTGGCCGCTGAGACCCTCGCGCCCCCGCCCCC
TCTTTTCCCCTTGATGGAATTGACCATAACAATTTTCATCTCCTTTCAGGGGATCAAAGGACGGAGTG
GGGGACAGAGACTCAGATGAGGACAGAGTGGTTCCAATGTGTTCAATAGATTTAGGAGCAGAAATGC
AAGGGGCTGCATGACCTACCAGGACAGAACTTCCCAATTACAGGGTACTCACAGCCGATTGGTGA
CTCACTCAATGTGTCATTTCCGGCTGTGTGTGAGCAGTGGACACGTGAGGGGGGGTGGTGAGA
GAGACAGGCAGCTCGATTCAACTACCTTAGATAATATTTCTGAAAACCTACCAGCCAGAGGGTAGGGC
ACAAAGATGGATGTAATGCACCTTGGGAGGCCAAGGCGGGAGGATTGCTTGGCCAGGAGTTCAAGAC
CAGCCTGGGCAACATACCAAGACCCCGTCTCTTAAAAATATATATATTTAAATATACTTAAATATA
TATTTCTAATATCTTAAATATATATATATATTTTAAAGACCAATTTATGGGAGAATTGCACACAGATG
TGAAATGAATGTAATCTAATAGAAGCCTAATCAGCCCACCATGTTCTCCACTGAAAAATCCTCTTTCTT
TGGGGTTTTCTTCTTTCTTTTGGATTTGCACTGGACGGTGACGTCAGCCATGTACAGGATCCACA
GGGGTGGTGTCAAATGCTATTGAAATTGTGTTGAATTGTATGCTTTTCACTTTTGATAAATAAACATG
TAAAAATGTTCAAAAAATAATAAAATAAATAACGAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_000602.5](#)

Summary: This gene encodes a member of the serine proteinase inhibitor (serpin) superfamily. This member is the principal inhibitor of tissue plasminogen activator (tPA) and urokinase (uPA), and hence is an inhibitor of fibrinolysis. The protein also functions as a component of innate antiviral immunity. Defects in this gene are the cause of plasminogen activator inhibitor-1 deficiency (PAI-1 deficiency), and high concentrations of the gene product are associated with thrombophilia. [provided by RefSeq, Aug 2020]

Locus ID: 5054

MW: 69.2