

Product datasheet for **SC216720**

FBXO33 (NM_203301) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FBXO33 (NM_203301) Human 3' UTR Clone
Symbol:	FBXO33
Synonyms:	BMND12; c14_5247; Fbx33
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_203301
Insert Size:	1862 bp



[View online »](#)

Insert Sequence: >SC216720 3'UTR clone of NM_203301
The sequence shown below is from the reference sequence of NM_203301. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAGATGCAAAGCTTCAGTGAAGACATTAGCTTTTTTTTAAATGTAGAATTCTGTGGTTACATATGCAA
GTAGGGTCTATTATGTTTTTTTTCAGTAGTGTGAATTAATCCTTTTGTGCTGTGTTTAAATCAGTATT
AGCTTTATAGAATTATATATGTATATTCTACTTCTTGATCAAAGAACGTAGTCGGGTATTGGTTTAGAA
GTTCAAAGTGACAATGTATAGGGCTTTACGGTTAATGGACTTGTACCAAACCTTAAGGATATACAGC
CGAAGATTGTCTGAGGTGTGGCTAACTTTATTTTCACTGAGTTACTCTGCCTTTTACGTTTTTA
TTCTTTGTGTGCAGAGTTCAGAGCTCAGGAGCCAAAATATTTTATACATATATAGATATATATCCAT
AGCCTGGTAGATTATATGCAATGCACTGCATCCATGCATTCTGATAGCATTTTCATTAATTTTGATTTC
AAATGGAAGATAAATAATATCCCAAATGAATTATCTGTAAACAGAAAAGCCAAGACTTTAACTTTTCATTA
CATATCTAATAGTTGATATCACCAGTTACCATTTTGAATTTTGTATAGTACTAGGTTAGAACATTGCCTT
AATCCTTTTAAAAAAATGCATTTACGTAAACAGAACTACTGAAATTGTTGGATTTTTTAACTATATCT
GACATAATTTTATTCATCAATTACATTACATTCATTTAGCCCAGATTTCAAATAGTGGGGGAAGAAA
GAAACTGTATTTAGAGTAAAATCTCCTAAAGGAAATAAAACACAGAGTTGTAATACACATGCTTGC
AAAAACATTAGTCGTGAAATCCCTAGCAACAAGTCACTGGATTTTTCTCTGTGCAGCAGCGGTGCAGCT
GCCAAAGAATAGACTTAATGAAGAAGTGGCCACATGCTGGCAGGGGCCGGCCACTCTGGCCAGCCAGA
TACTGCTAGATTGTAATATTTAAGGTGCAATTTGACCTGTGGTACACAGCTGTGCTGTGGCTCAGTCA
GCAACCTCAGAACTCTGAAAAACAAAAACAAAAACAAAAACAAAAACAAAAACAAAAACAAAAACATGCACCTGTTT
CACTGTGAATAGTGAATGTAAAGGAAAGAAAGGAAAAACCTGAAAGCTTGTCTATCACAGGTATGAGCT
GCTATGATACATGAAGAACATTCATGAAGTATGTTTTAAACCTTGTTATATCTGAGAGCTTTAAAA
GCCTATTAAGTGTTCAGGGCAACCGCGGTACAGAGCTGGTCTCTGTGAGACTTCCACCTGACCCAGT
TTTAAGTGGTACGAATGTTGTGCATTTAATGTTAAAGGACAGTCTGCAATAATAAGTAAGTAGCCAAC
CTGGGTGCCAGCAGTGTGAGACCTGGCTGCTCTATTGTAAGCTTTGGAACACAATTTATGCAACAG
ATGTCCAGATATGATTCTATTTATGAAAAAGTTTATATGTGTTTTACAAATGGTTTTACCATCTTAT
ATTAATAAGACCTTTTACAGGTGTGCACTGTTTTGTCTCCAGTGAGCACATACCATGCGGATTTTATA
TGTACATCAGTAGTGTGAATCCACTGGCACAGTGTGTGTAATGCCAGATGTGGTGAATTTTATCTTA
TAAATGTGATCAGATAAAATAACTCCTGACAGAACTGTAAAGGAAACCAGCTGAATGTTTGACCTGATG
ACTGATGTTGTATGGTTTATGTTAAATGTATATCTTTTAAATCAATGAATAAAGCATTAAAAATGATCA
TTGTAATATATTTTATTGTATATGAAGCATGACTACTTTTGTAGTTAAATAAACTGCAATAATCCTTA
AGCGGACCGACTTACGCGTAAGCGGCCGCGGCATCTAGATTGCAAGAAAAAGACCGACCAAGCGACGCC
CAACCTGCCATCACGAGATTTTCGATTCCACCGCCGC
  
```

Restriction Sites: Sgfl-RsrII

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_203301.4](#)

Summary: This locus represents an member of the F-box gene family. The encoded protein contains an F-box motif and a domain that might form a structure similar to a leucine-rich repeat found in placental RNase inhibitor. This locus may be associated with copy number variation of UGT2B17 (GeneID 7367), which has been associated with susceptibility to osteoporosis. [provided by RefSeq, Sep 2010]

Locus ID: 254170

MW: 72.3