

## Product datasheet for **SC216726**

### **CANT1 (NM\_001159772) Human 3' UTR Clone**

#### **Product data:**

Product Type:	3' UTR Clones
Product Name:	CANT1 (NM_001159772) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CANT1
Synonyms:	DBQD; DBQD1; EDM7; SCAN-1; SCAN1; SHAPY
ACCN:	NM_001159772
Insert Size:	1864 bp



[View online »](#)



**Summary:** This protein encoded by this gene belongs to the apyrase family. It functions as a calcium-dependent nucleotidase with a preference for UDP. Mutations in this gene are associated with Desbuquois dysplasia with hand anomalies. Alternatively spliced transcript variants have been noted for this gene.[provided by RefSeq, Mar 2010]

**Locus ID:** 124583

**MW:** 68.9