

Product datasheet for **SC216749**

SLC25A35 (NM_201520) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLC25A35 (NM_201520) Human 3' UTR Clone
Symbol:	SLC25A35
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_201520
Insert Size:	1883 bp



[View online »](#)

Insert Sequence: >SC216749 3'UTR clone of NM_201520
The sequence shown below is from the reference sequence of NM_201520. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGCCCGATTCTCCTGGGAACCCAGTCAAGATACCCGCCCGTTCCCTCCTTACGCGGGGCTGATTG
AAGGTAAGCAAGAGATCAGTCTGGTACTGGGGCAGCCTCAGCAAGGCCTGATGAAGTGTCACGTCCTTT
GCTACCTAAGGGGATCAGAGTCAGGATGCTTACCTGGATCTCTGCGTCCCGCTACCGCCCTTCCAGGA
CATTACACACTCTCGGGCCCTCACCTGCTGGTTTTCTTAGCTATCTGCTGCTTGCAGAGAGGACCCA
GGCTTCTTGACAGCGGCCCTCAGGGCCAGGTTCTCAAAGTGGTACCTGGAGTAGGGAAGCGGGTCAAGAA
AGCCCTAGCCCCCTGCACGCCACCAACATCCTCAAAGTGGTACCTGGAGTAGGGAAGCGGGTCAAGAA
GGGACATCTCTGGTGCCTCCCGGGAGCATCTTAAAGTTGCAGGCCTGGCTCCCGGCCCTGCGTGGT
CTACCCGCCCGGATTAACCGTCGCGCCGCGGGCCCGCCACTGGTCCGACCCCGCCCTGCCATTGGCC
AGCCGGGGGCCATTCCCTCACCGCGCAGCGCTTCGCCCCGTACGTGGGCTGCAGCTCGAGAAGTTCC
ACTATCAGGCTCTGGTCCGTACGGGATGGCAGAAAATCTTGTATTGTCCGGGACCGGTGCGAGGTGCG
CTGGAGTAGGGCAGGCAGACCTGTTGACCCCGGTCTCTTCCACAACCCACCGAGTACGCGCCCCC
TGGGCGCCCCGGCCTTCTACCTTACGTCAATGGCCCCCATGGGGAGGATGGCGGAAAAGCGCCCCCG
AACAGCGGGCAGTCTCTCGTGGGCTCCATGGGTCTGAGGCTGGGTATAGGGCATTAGAAAGGTTAAAA
GAAGTCCACATTGGGTACGTGGCATCGGCTCCCTGGGTTTCGAGCTGCCACCGCCACCGGGTCTGG
CTTAGCTCCCGGATCCTTAAAGATTGGCTCCGCCCTGCGCTCTCGTTGGCCAATCGGCGCAAGACCC
TGGGGCGGGGGTCAGCGGCTTTGCTCCGGAAGCCGCGCTCTTAGGAGTGCCTGCGGGGACGGTGGG
GTCTGCCTTGGGGTTGGACTTAGGGTTAGGGTCTATGTCGCCCCCGGTGGTACAGAGCATGCTCCCGCC
AGCCTCCTCTTGATTGATTGATTGATTGATTGATTGACAGGTCTCGTGTGTCGCCAGGCTGTAATGC
AGTGGCGCATCACCGCTCACTGCACTCTCACCTCCAGGACTCAAGTGATACTCTGCCTCAGCCTCT
TGAGTGGCTCTGCTCCCGCCACCACGCCCGGCTAAATTTTAACTTGTAGAAATGGGGCGGGTG
GGGAGTGTCTCCAATCTTGCCAGGCTGGCTCAAAGTCTGGACTCAAGCGATCTCCCGCTCGGC
TCCCTAAGTGCTGGAATTACAGGCGAGAGCCACAGCGCTGGCTCCTCGTCTTTCAAGTGGTCACA
CCTTCAAAGAAGTGCCTGAGACCTCTGAGGGGAGGTGCTCGTCAAAGGGGCCACGTCTAGGTA
GGGTCTCTTGAAGTCAAGAAAGCATAGTTTACGGTGTCTGGGCTCTGAAAGAAGATGGGCTGCCAC
CCAGACCGCTAATTCATAGATCTGCTGCTGGGAATTAGGGGAACAATGCGTCCAGGTTAGGATTCAGGT
CTCACATCTTCTGAAAGTGAAGCCAACTGGGTATTAGAAATCTAATTCCTCTCAAATATTTGCCTT
GGGAATGCTCCCGGGCAGGTGCAACCCACCCGAGCTCCACACCTGGACAGAAGCACAGATTAAC
CCTCTTTGTTTGTCTTTA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_201520.3](#)

Summary: SLC25A35 belongs to the SLC25 family of mitochondrial carrier proteins (Haitina et al., 2006 [PubMed 16949250]).[supplied by OMIM, Mar 2008]

Locus ID: 399512

MW: 66.7