

Product datasheet for **SC216817**

Versican (VCAN) (NM_001164098) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Versican (VCAN) (NM_001164098) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	VCAN
Synonyms:	CSPG2; ERVR; GHAP; PG-M; WGN; WGN1
ACCN:	NM_001164098
Insert Size:	1898 bp



[View online »](#)

Insert Sequence: >SC216817 3'UTR clone of NM_001164098
 The sequence shown below is from the reference sequence of NM_001164098. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGCCGGAGGTGGCAGGAGTCGAGGCGCTGATCCCTAAAATGGCGAACATGTGTTTTTCATCATTTTCAGCC
AAAGTCCCTAACTTCTGTGCCTTTCCTATCACCTCGAGAAGTAATTATCAGTTGGTTTGGATTTTGGGA
CCACC GTTCAGTCATTTTGGGTTGCCGTCTCCAAAACATTTTAAATGAAAGTATTGGCATTCAAAAA
GACAGCAGACAAAATGAAAGAAAATGAGAGCAGAAAAGTAAGCATTTCAGCCTATCTAATTTCTTTAGT
TTTCTATTTGCCTCCAGTGCAGTCCATTTCTAATGTATACCAGCCTACTGTACTATTTAAATGCTCA
ATTTTCAGCACCGATGGCCATGTAATAAGATGATTTAATGTTGATTTTAACTCTGTATATAAAATAAAA
AGTCACAATGAGTTTGGGCATATTTAATGATGATTATGGAGCCTTAGAGGTCTTTAATCATTGGTTCGG
CTGCTTTTATGTAGTTTAGGCTGGAATGTTTTCACTTGCTCTTTGACTGTCAGCAAGACTGAAGATGG
CTTTTCCTGGACAGCTAGAAAACACAAAATCTTGTAGGTATTGCACCTATCTCAGCCATAGGTGCAAGT
TTGCTTCTACATGATGCTAAAGGCTGCGAATGGGATCCTGATGGAAC TAAGGACTCCAATGTGCAACTC
TTCTTTGCTGCATTCTTTTTCTTCACTTACAAGAAAGGCCTGAATGGAGGACTTTTCTGTAAACCAGGA
ACATTTTTTAGGGGTCAAAGTGCTAATAATTAACCAACCAGGTCTACTTTTTAATGGCTTTCATAACA
CTAACTCATAAGGTTACCGATCAATGCATTTTCATACGGATATAGACCTAGGGCTCTGGAGGGTGGGGGA
TTGTTAAAACACATGCAAAAAAAAAAAAAAAAAAAAAAAAAAAGAAATTTTGTATATAAACCATTTAAT
CTTTTATAAAGTTTTGAATGTTTCATGTATGAATGCTGCAGCTGTGAAGCATACATAAATAAATGAAGTA
AGCCATCTGATTTAATTTATTGGATGTTATTTTCCCTAAGACCTGAAAATGAACATAGTATGCTAGTT
ATTTTTCAGTGTTAGCCTTTTACTTTCCCTCACACAATTTGGAATCATATAATATAGGTACTTTGTCCCT
GATTAATAATGTGACGGATAGAATGCATCAAGTGTATTATGAAAAGAGTGGAAAAGTATATAGCCTT
TTAGCAAAAGGTGTTTGCCCATTTCTAAGAAATGAGCGAATATATAGAAATAGTGTGGCATTCTTCCCT
GTTAGGTGGAGTGTATGTGTTGACATTTCTCCCATCTCTTCCCACTCTGTTTTCTCCCATTTTGA
ATAAAGTGACTGCTGAAGATGACTTTGAATCCTTATCCACTTAATTTAATGTTTAAAGAAAAACCTGTA
ATGGAAAGTAAGACTCCTTCCCTAATTTTCAGTTTAGAGCAACTTGAAGAAGAGTAGACAAAAATAAAA
TGCACATAGAAAAAGAGAAAAAGGGCACAAAGGGATTGGCCCAATATTGATCTTTTTTTATAAAACCT
CCTTTGGCTTAGAAGGAATGACTCTAGCTACAATAATACACAGTATGTTAAGCAGGTTCCCTTGGTTG
TTGCATTAATGTAAATCCACCTTAGGTATTTAGAGCACAGAACAACACTGTGTTGATCTAGTAGGTT
TCTATTTTTCTTTCTTTTACAATGCACATAACTTTCTGTATTTATATCATAACGTGTATAGTGT
AAAATGTGAATGACTTTTTTTGTGAATGAAAATCTAAAATCTTTGTAAC TTTTATATCTGCTTTTGT
TCACCAAAGAAACCTAAAATCCTTCTTTACTACA
ACGCGT AAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001164098.2](#)

Summary: This gene is a member of the aggrecan/versican proteoglycan family. The protein encoded is a large chondroitin sulfate proteoglycan and is a major component of the extracellular matrix. This protein is involved in cell adhesion, proliferation, proliferation, migration and angiogenesis and plays a central role in tissue morphogenesis and maintenance. Mutations in this gene are the cause of Wagner syndrome type 1. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2009]

Locus ID: 1462

MW: 74.3