

## Product datasheet for **SC216835**

### LANCL3 (NM\_198511) Human 3' UTR Clone

#### Product data:

|                           |                                       |
|---------------------------|---------------------------------------|
| Product Type:             | 3' UTR Clones                         |
| Product Name:             | LANCL3 (NM_198511) Human 3' UTR Clone |
| Symbol:                   | LANCL3                                |
| Mammalian Cell Selection: | Neomycin                              |
| Vector:                   | pMirTarget (PS100062)                 |
| ACCN:                     | NM_198511                             |
| Insert Size:              | 2000 bp                               |



[View online »](#)

**Insert Sequence:**

>SC216835 3'UTR clone of NM\_198511

The sequence shown below is from the reference sequence of NM\_198511. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

[illegible]

**Restriction Sites:**

Sgfl-MluI

### OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

### Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

## RefSeq:

NM 198511.3

**Locus ID:**

347404

MW: 76.6