

Product datasheet for **SC216848**

B3GALNT1 (NM_033168) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	B3GALNT1 (NM_033168) Human 3' UTR Clone
Symbol:	B3GALNT1
Synonyms:	B3GALT3; beta3Gal-T3; galT3; Gb4Cer; GLCT3; GLOB; P; P1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_033168
Insert Size:	1906 bp



[View online »](#)

Insert Sequence: >SC216848 3'UTR clone of NM_033168
The sequence shown below is from the reference sequence of NM_033168. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGCATCGCC
ATGCTAAGGAACACCACATGCCATTATTAACTTTCACATTCTACAAAAAGCCTAGAAGGACAGGATACTT
TGTGGAAGTGTAAATAAAGTAGGTACTGTGGAATAATCATGGGAGGTCAGTGTGCTGGCTTACACT
GAACTGAACTCATGAAAAACCCAGACTGGAGACTGGAGGGTTACACTTGTGATTTATTAGTCAGGCC
TTCAAAGATGATATGTGGAGGAATTAATATAAAGGAATTGGAGGTTTTGCTAAAGAAATTAATAGGA
CCAAACAATTTGGACATGTCATTCTGTAGACTAGAATTTCTTAAAGGGTGTACTGAGTTATAAGCTC
ACTAGGCTGTAAAAACAAAACAATGTAGAGTTTTATTTATTGAACAATGTAGTCACTTGAAGTTTTGT
GTATATCTTATGTGGATTACCAATTTAAAAATATATGTAGTTCTGTGTCAAAAAACTTCTCACTGAAG
TTATACTGAACAAAATTTTACCTGTTTTTGGTCATTTATAAAGTACTTCAAGATGTTGCAGTATTTAC
AGTTATTATTATTAAAAATTACTTCAACTTTGTGTTTTTAAATGTTTTGACGATTCAATACAAGATAA
AAAGGATAGTGAATCATTCTTTACATGCAACATTTTCCAGTTACTTAACTGATCAGTTTATTATTGAT
ACATCACTCCATTAATGTAAAGTCATAGGTCATTATTGCATATCAGTAATCTCTTGGACTTTGTTAAAT
ATTTTACTGTGGTAATATAGAGAAGAAATTAAGCAAGAAAATCTGAAGTATTGTCTTGTTTTTAAAAA
TACAGTTCCTAGTGTTTTAGAAGTCACCTAATTTGTCTCATTTTTCCACCTGGAATTAGGAATAATG
TAGAATGCAAGGCAGTAATTTCTTTTGGAAAGGACTCTGAAGGCAGAAAAGAAGGGAGAGAACCTCAT
GGGCAGAAATATTATAAAAAAGAGTGTATATCCAGCATTTGAATTGGAAGAGAGAGTGAAGATCCAA
GTTGCATTATTAATCTGCCCTGTGTTTTTCTTTTAAACAATCAGTTTGAGCTGCTGCTGTTATGAGTT
TCTCATCAAGATGAAAGCCCTAATATGTAAAGTCAAAATCCGATTTAAATTTTGTGCTTTTATAGAAAGA
AATTTCTTCATAGACGTGGTGATATATCATTTGTTGGACCTGCTAATAGTAGGTCAAAGGGAGCACTC
CTTGCCCCCTGTTCTGGGTTTATGCAGTTTTCTTTTATAGAGTTTATATAGGGCAAGTGTTCTTTTTC
TCTGAATTACAGGATGGAAGGTCATATCCTTTGTCAGGAAATATAAACTTGAAGTATGTAGTCAG
CTCTTGTAATACTCATATTTATGATTGTCCTATATGAAAAACAACTTCAGTTAAACTATAATGTGTGA
TTCTGTATAACAAGGTGATGTCTGTTTCCAGGGCTCAGACCTAATCCAGTTATAATAAAATCAATTAA
ATGAAATATTCTATAGAATCGATCTATGCCCTTGTTAATCTCCATCCATATAGGAGTCACGTTCTTTAA
GACAGATGGTGGTAGTTATTTTGTGGCATGGTTAGATTTGACTGGTTTTGCAGAAGATTACAGTTATG
TACTGCATAATGACATATACAATAGTGGTCCCATAAAATTATAATGGAGCAGAAAATCTATTGCCTCAT
GATGTTTTAGCCGTTTTAATGTCATAGCCTAGTGCATTACTCACGTGTCTGTGGAGATGTGGTGTAAA
CAAACCTACTACACTGCCAGTTGTATAAAAGTACAGCACATTCAGTTATGTACAGTATGTAATACTGAT
AATGACAATAAATGACACCGGTTTGTGTATTTACTTTTTATAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_033168.3](#)

Summary:

This gene is a member of the beta-1,3-galactosyltransferase (beta3GalT) gene family. This family encodes type II membrane-bound glycoproteins with diverse enzymatic functions using different donor substrates (UDP-galactose and UDP-N-acetylglucosamine) and different acceptor sugars (N-acetylglucosamine, galactose, N-acetylgalactosamine). The beta3GalT genes are distantly related to the Drosophila Brainiac gene and have the protein coding sequence contained in a single exon. The beta3GalT proteins also contain conserved sequences not found in the beta4GalT or alpha3GalT proteins. The carbohydrate chains synthesized by these enzymes are designated as type 1, whereas beta4GalT enzymes synthesize type 2 carbohydrate chains. The ratio of type 1:type 2 chains changes during embryogenesis. By sequence similarity, the beta3GalT genes fall into at least two groups: beta3GalT4 and 4 other beta3GalT genes (beta3GalT1-3, beta3GalT5). The encoded protein of this gene does not use N-acetylglucosamine as an acceptor sugar at all. [provided by RefSeq, Mar 2017]

Locus ID:

8706

MW:

75.5