

Product datasheet for **SC216881**

ZHX1 (NM_001017926) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZHX1 (NM_001017926) Human 3' UTR Clone
Symbol:	ZHX1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001017926
Insert Size:	1898 bp



[View online »](#)

Insert Sequence:	>SC216881 3'UTR clone of NM_001017926 The sequence shown below is from the reference sequence of NM_001017926. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC AAACGGAAGCTGTCTAAATCAGATGACGAAATCTGCCTAAACGTTGAAGGAGAATCAATCTTCAAC TCAAGATGTCTGATTTACTGTGAATTTGCCAATCTTTGATGACATTGAAAACGTTTTGGGGCATACAC ACTCAAAAAGCAGGATCCAATACCCAAAAGAAATGGAACCTTAATGTTGTGCCAAAGTTAACTACTGCA GTTGGTGAAGTTCTGCAATGTAATAGAACACTAATTAACAACTTGTAAAAATGCAATTTAAAT TTAATACAGTACATTTTCTTCTAATATGATGGAGACATTCTGAATCTTAGACTTTCTGAGGGGTTT AATGACCACTAGAGCTTGTCTCATATTCAGTCCAGTTTAATACTGTATGTCTAGTAAGATGGGCTATA TATTGCCTCTATTCTTTGAGATGTGATTAAGCTTAGAAGTTTACCAAAAAATTTACAGTAAAAATGT TAGAAATGGGTAAAAAAGGCTGATTTTACTATGTCTTATTTATCAGACCT GCACTAAGATTAAAGTTGTGCATGGGCTTATACAATTTCAAATGTAAAAACAGCCTATTATATATGT TGAAACAATAAGATTTTATGGTTTCACTGGCCCTTAATAATTGGATCCTTCATTAACACATCCATGAA TACCACTGTTTCTACTGTTACTTGTCCATTAGAAGCTGATTTGCTTTTATTAACCTGCAGGAAATTA ATTCTGAAGAAGGAAAGAACTGACATGGTAATTTGGATAGGGTGAAGCTTTGAAATAAATACTGCAC CAAAAATGTGAAAAGAATCACGTAGTAAATGTATACTTAGTGTTTTCTTACAATCGTGGAAATCAGAA ATGTTTATACAGGTGTCAGATCAGAATTTTAAAGATTATCTTTCAGTGAAATTTTATAACAGTTTCAG TGAATAGAGAAAAATTTGAACTGATCATAATGAGTTAAAGTTAGCTTTTCACTTATGAATCATGGG AGTGTTAACTGTTAACTCTTGCGTATCAGGAGACTGAAAGTTCACTAATTTACAACTTTGTAAAG ACTGGTGCTGTTTCAATCATATTGGGAAGAACCCTTTCTTTGTGACCATAGGACCATTTTCAAAATGC GCTTTCTCGTAGAAGAAATGTTCAAGTTTAATTAGCTTTGGATTTTAACTGTTTTCTGAATGACCAG ATGTAGGCTTAGCCAGTTTATTTTCTCTAATTTAATGGATAATAAATATTAAGCTCTTGTTAAATGTT TTTCAGCCATTGAAAATGTTATGAAAAAGCCTTGAATATGCATGCTGCTTTTATTTTATAAACTTTT ATTTTTAACTATAGCTTTATTAGTAATCCAAAATGCTGCAATGTAGCTGATTTCTCTACTCAGACAG CTGGCTTTGTGGTGGCTTTTTCATACTACTGTTAACTTTTTAAAAAGAAGCAATGTAAAGACCGT AATGCACTAAGCTTTTTTTTTTCTTTTACACGTTTAAAACTTCTTAAAGCACTCTGTAATATAATGA TTAAATTGCTCCCTTTTATAACATTGCTAATGTGGTCTGAGTTTGTGTTGACAATAATTTAATGTTC TCAAAATTTGAAATGTCAAGTATTAGAAATGAAGCCTTTGAATATAGTTGTTCCCTTGCAATATGTTA ACTCTGTGTAGCTGTATATACTATTAATTTCCCTCATCTCTTTATAGGCTGTAAAACTCTTCTCT GACGTTTGCTAGCTTTTAAATATATTTATTTGACAGAGCAAAAGAACCCTCTTCTGGTTGATTGTGACC TTGAATACAAAATAAAGTAGTGATTAGTAGTTTA ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001017926.3</u>

Summary:

The members of the zinc fingers and homeoboxes gene family are nuclear homodimeric transcriptional repressors that interact with the A subunit of nuclear factor-Y (NF-YA) and contain two C2H2-type zinc fingers and five homeobox DNA-binding domains. This gene encodes member 1 of this gene family. In addition to forming homodimers, this protein heterodimerizes with members 2 and 3 of the zinc fingers and homeoboxes family. Alternative splicing results in multiple transcript variants. Read-through transcription also exists between this gene and the downstream chromosome 8 open reading frame 76 (C8orf76) gene. [provided by RefSeq, Feb 2011]

Locus ID:

11244

MW:

74.7