

Product datasheet for **SC217250**

RAB37 (NM_001163989) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RAB37 (NM_001163989) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	RAB37
ACCN:	NM_001163989
Insert Size:	1954 bp



[View online »](#)

Insert Sequence: >SC217250 3'UTR clone of NM_001163989
 The sequence shown below is from the reference sequence of NM_001163989. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAGCGCTCCAGCTGCTGCTCCTTCATGTGAATCCAGGGGCAGAGAGGAGGCTCTGGAGGCACACAGG
ATGCAGCCTTCCCCCTCCAGGCCCTGGCTTATTCCAAGAGGCTGAGCCAATGGGGAGAAAGATGGAGGA
CTCACTGCACAGCCGCTTCTAGCAGGGAGCTATACTCCAACCTCTACTTGAGTTCTGCGGTCTCCCC
GCATCCACAGGGAGGGTAAAACACTTAGCTTTTATTTTAAAGTACATAATTTAATACAAAAAAGGCG
CCTGGATCCCCAAAAACCGAGGCTGGGAGCTAGTGGCCCTTTTGTCTTCTAGGACTTGGGGGGCCGGC
CCTCCCTCCTAAGCATAACAAAGGTGGTGTGCTCCAGCTCAGCCCCAGGGGACACAGATGCACTTTGG
GGGTGAGGGCAGGTAATGACTCCATCGCACCTCAGTTCAGCTGGACAGAGGCTCAGGTGACCCAGCC
TTCAGTGTCTCCGCTCTCCAGGAGCTTATCTTCCGCCATCTCCAAATAAGTGGGCCCTTGCTGTG
GAGGAAGACCAAAGCCTCAGGGAAGATAAGAGATATGGAGATGGGAGGGGAGGACAAGGGGCAGAGAG
TAGGGTCTAGCTGGCTATCTCTGGCCTTACTAACACCCCTGGAGGCATGCCCTTTTCTCCAGCACA
CAAGCACATTGGGGCACCTGAAATATTGGTTCAGGCTCCTGTTCTCTGGACTTCAGATCCTGGGGGA
GCCCTCCCCCCTGAATCCCTGGCTTAGCTACCTTCTGCCTGTGCACCTAAAAACCTCAGGTGAGA
ACTAGGAAAAGAGTTTGTTTTTATTTTTTGAATGGAGTCTCGTTCTGTGCGCCAGGCTGAGGTGCA
GTAGTGAATCTCCGCTCACTACAACCTCCACTCCCTGGGGCTCAAGCGATCCTCCACCTCAGCCGCC
GAAGTACTGGGACTATAGGTGTGTACCATCACACCTGGCTAATTTTTGTATTTTTGTAGACACAGGG
TTTCGCCATGTTGCCAGGCTGGTCTTGAATCCTGAGCTCAAGCAACCTGCCGGCTCGGCCCTCCAA
AGTACTGGGATTACACGCAGAAGGCACCATGCCAGGCTAGATGTGCTTATCCCAATCCTTTGGCAGG
CATGCAGCTCCACAGGCGATTTCTCAAGCAGCTGAAGTGTTTAGCCCTCCTGGGTTAAGAGCCAGATA
AGGAGAAATCCCTTCTAGGTTTGAATGTGTTGAAAAAAGAGAAATCCCTGGCTCCTGGAGCT
GGTGGGAGACAAGATTAAGCAAACCTCCCTGACATGTATCCCTTTGACCCCAAGCTCTGCCTCCTCCC
TGACCACCATGCCCTTCTTTAACTTCTCAAACAGATACCAGGGCCTAAACTGCTTTACCTCCCTC
CTACTGAGTCAGGTTAGGTGGTGGGAGGTCACCCATTTCCGAGTTAAACCAATGCAATATGAGTAAAC
AAAGTCATGTGGGTATGTCTGGGTAGAGAGAGGGGTAGCAAGTTCATGTGCTCCTCTTGGTCACATAT
CTCCAAAGCTCCGATCCCTGCCATGGGAAGTGGACAGGAAACATGAGGTCATGACTGCAGGCATCTT
TACTGCAGCTCTGCCGGCTGGAGGGGGAGAGGGGGAGGAAGAAGTATGCGCTGCACATTTCTGAGGCT
ACTGCATTTGCTTCAAGGCAGAAATCTTGTCTGAGCAGTCAGCGGCTCCAGTTTGGGGCCGATAAGG
AAGTTCTCCGTGGCTCCCTCAGGCAGAGCAGGGAGGAGGCTGACATTGCCAGTCTCTTCTGGGGCCCA
AGGCAGGTTGCAGGAGATCCAATCCCATAGACAGCTCTGGGCCTCTTGCAATTTGAGTTTTTCAGAATTA
AACTGCAGTATTTTGGAAAGCA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001163989.2](#)

Summary: Rab proteins are low molecular mass GTPases that are critical regulators of vesicle trafficking. For additional background information on Rab proteins, see MIM 179508.[supplied by OMIM, Apr 2006]

Locus ID: 326624

MW: 71.4