

Product datasheet for **SC217298**

DCAF12 (NM_015397) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	DCAF12 (NM_015397) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	DCAF12
Synonyms:	CT102; KIAA1892; TCC52; WDR40A
ACCN:	NM_015397
Insert Size:	1993 bp



[View online »](#)

Summary: This gene encodes a WD repeat-containing protein that interacts with the COP9 signalosome, a macromolecular complex that interacts with cullin-RING E3 ligases and regulates their activity by hydrolyzing cullin-Nedd8 conjugates. [provided by RefSeq, Jul 2009]

Locus ID: 25853

MW: 74.6