

Product datasheet for **SC300028**

Cytokeratin 9 (KRT9) (NM_000226) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cytokeratin 9 (KRT9) (NM_000226) Human Untagged Clone
Tag:	Tag Free
Symbol:	KRT9
Synonyms:	CK-9; EPPK; K9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC300028 representing NM_000226.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAGCTGCAGACAGTTCTCCTCGTCTACTTGAGCCGACGCGGGGGTGGCGGGGGCGGCCTGGGC
AGCGGGGGCAGCATAAGGTCTTCTACAGCCGTTACAGTCTCAGGGGGCGGTGGAGGAGGGGGCCGA
TTCAGCTCTTCTAGTGGCTATGGTGGGGGAAGCTCTCGTGTCTGTGGAGGGGAGGCGGTGGCAGTTTT
GGCTACAGCTACGGCGGAGGATCTGGGGTGGTTTTAGTGCCAGTAGTTAGGCGGTGGCTTTGGGGT
GGTTCCAGAGTTTTGGTGGTCTTCTGGAGGAGGCTATAGTAGTTCTGGGGTTTTGGAGGTGGCTTT
GGTGGTGGTCTGGAGGTGGCTTTGGTGGTGGCTATGGGAGTGGTTTTGGGGGTTTTGGGGCTTTGGA
GGTGGTGTGGAGGAGGTGATGGTGGTATTCTGACTGCTAATGAGAAGAGCACCATGCAGAACTCAAT
TCTCGGCTGGCCTCTTACTTGGATAAAGTGCAGGCTCTAGAGGAGCCAAACAACGACCTGGAGAATAAG
ATCCAGGATTGGTACGACAAGAAGGACCTGCTGCTATCCAGAAGAACTACTCCCCTTATTATAACT
ATTGATGATCTCAAGGACCAGATTGTGGACCTGACAGTGGCAACAACAAAACCTCCTGGACATTGAC
AACACTCGCATGACACTGGATGACTTCAGGATAAAGTTTGAGATGGAGCAAAAACCTGCGGCAAGGAGTG
GATGCTGACATCAATGGCCTGCGGCAGGTGCTGGACAATCTGACCATGGAGAAGTCTGACCTGGAGATG
CAGTATGAGACTCTGAGGAGGAGCTGATGGCCCTCAAGAAGAATCATAAGGAGGAGATGAGTCAGCTG
ACTGGGCAGAACAGTGGAGATGTCAATGTGGAGATAAACGTTGCTCCTGGCAAAGATCTACCAAGACC
CTCAATGACATGCGTCAGGAGTATGAGCAGCTCATTGCTAAGAACAGAAAGGACATCGAGAATCAATAT
GAGACTCAGATAACCCAGATCGAGCATGAGGTATCCAGTAGTGGTCCAGGAGTGCAGTCCAGTGCCAAG
GAGGTGACCCAGCTCCGGCACGGTGTCCAGGAGTTGGAGATTGAGCTGCAGTCTCAGCTCAGCAAGAAA
GCAGCTCTGGAGAAGAGCTTGAAGACACGAAGAACCGTACTGTGCCAGCTGCAGATGATCCAGGAG
CAGATCAGTAACTTGGAGGCCAGATCACTGACGTCCGGCAAGAGATCGAGTGCCAGAATCAGGAATAC
AGCCTTCTGCTCAGCATTAAAGATGCGGCTGGAGAAGGAAATCGAGACCTACCACAACCTCCTTGAGGGA
GGCCAGGAAGACTTTGAATCCTCCGGAGCTGGAAAAATTGGCCTTGAGAGTTCGAGGAGGAAGTGGAGGC
AGTTATGGAAGAGGATCCAGGGGAGGAAGTGGAGGCAGCTATGGTGGAGGAGGAAGTGGAGGTGGCTAT
GGTGGAGGAAGTGGTCCAGGGGAGGAAGTGGAGGCAGCTACGGTGGAGGAAGTGGTTCTGGAGGAGGT
AGTGGAGGTGGCTATGGTGGAGGAAGTGGAGGTGGCCATAGCGGAGGAAGTGGAGGTGGTCATAGTGA
GGAAGTGGGGCAACTATGGAGGAGGAAGTGGCTCTGGAGGAGGAAGTGGGGTGGCTATGGTGGAGGA
AGTGGGTCCAGGGGAGGAAGTGGAGGCAGCCATGGTGGAGGAAGTGGTTTTGGAGGTGAAAGTGGAGGC
AGCTACGGAGGCGGTGAAGAAGCGAGTGAAGTGGTGGCGCTACGGAGGAGGAAGCGGAAAATCATCC
CATTCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_000226

Insert Size: 1872 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_000226.3](#)

RefSeq Size: 2290 bp

RefSeq ORF: 1872 bp

Locus ID: 3857

UniProt ID: [P35527](#)

Cytogenetics: 17q21.2

MW: 62.1 kDa

Gene Summary: This gene encodes the type I keratin 9, an intermediate filament chain expressed only in the terminally differentiated epidermis of palms and soles. Mutations in this gene cause epidermolytic palmoplantar keratoderma. [provided by RefSeq, Jul 2008]