

Product datasheet for **SC300042**

Plakophilin 1 (PKP1) (NM_000299) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Plakophilin 1 (PKP1) (NM_000299) Human Untagged Clone
Tag:	Tag Free
Symbol:	Plakophilin 1
Synonyms:	B6P; EDSFS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC300042 representing NM_000299.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAACCACTCGCCGCTCAAGACCGCCTTGGCGTACGAATGCTTCCAGGACCAGGACAACCTCCACGTTG
GCTTTGCCGTCGGACCAAAGATGAAAACAGGCACGTCTGGCAGGCAGCGGTGCAGGAGCAGGTGATG
ATGACCGTCAAGCGGCAGAAGTCCAAGTCTTCCAGTCGTCCACCCTGAGCCACTCCAATCGAGGTTC
ATGTATGATGGCTTGGCTGACAATTACAACATATGGGACCACCAGCAGGAGCAGCTACTACTCCAAGTTC
CAGGCAGGGAATGGCTCATGGGGATATCCGATCTACAATGGAACCCCTCAAGCGGGAGCCTGACAACAGG
CGCTTCAGCTCCTACAGCCAGATGGAGAACTGGAGCCGGCACTACCCCGGGGACGCTGTAAACACCACC
GGCGCAGGCAGCGACATCTGCTTCATGCAGAAAATCAAGGCGAGCCGAGTGAGCCCGACCTCTACTGT
GACCCACGGGGCACCTGCGCAAGGGCACGCTGGGCAGCAAGGGCCAGAAGACCACCCAGAACCCTAC
AGCTTTTACAGCACCTGCAGTGGTCAGAAGGCCATAAAGAAGTGCCTGTGCGCCCGCCCTTTGTGCC
TCCAAGCAGGACCCTGTGTATATCCCGCCATCTCCTGCAACAAGGACCTGTCTTTGGCCACTCTAGG
GCCAGCTCCAAGATCTGCAGTGAGGACATCGAGTGCAGTGGGCTGACCATCCCAAGGCTGTGCAGTAC
CTGAGCTCCCAGGATGAGAAGTACCAGGCCATTGGGGCTATTACATCCAGCATACCTGTTCCAGGAT
GAATCTGCCAAGCAACAGGTCTATCAGCTGGGAGGCATCTGCAAGCTGGTGGACCTCTCCGAGCCCC
AACCAGAACGTCAGCAGGCCCGCGCAGGGGCCCTGCGCAACCTGGTGTTCAGGAGCACCACCAACAAG
CTGGAGACCCGGAGGCAGAATGGGATCCGCGAGGCAGTCAGCCTCCTGAGGAGAACCAGGAAACGCCGAG
ATCCAGAAGCAGCTGACTGGGCTGCTCTGGAACCTGTCTTCCACTGACGAGCTGAAGGAGAACTCATT
GCCGACGCCCTGCCTGTCTGGCCGACCGCTCATCATTCCCTTCTCTGGCTGGTGGCATGGCAATAGC
AACATGTCCCAGGAAAGTGGTGGACCTGAGGTCTTCTTCAATGCCACAGGCTGCTTGGAAGAGAGACTG
GGCATGCGGGAGCTTCTGGCTCTTGTTCGCAAAGGGCCACTAGTAGCAGGGTGAACCTGAGCTCGGCC
GATGCAGGCCGCCAGACCATGCGTAACTACTCAGGGCTCATTGATTCCCTCATGGCCTATGTCCAGAAC
TGTGTAGCGGCCAGCCGCTGTGACGACAAGTCTGTGAAAACCTGCATGTGTGTTCTGCACAACCTCTCC
TACCGCTGGACGCCGAGGTGCCACCCGCTACCGCCAGCTGGAGTATAACGCCCGCAACGCCTACACC
GAGAAGTCTCCACTGGCTGCTTCCAGCAACAAGAGCGACAAGATGATGAACAACAACCTATGACTGCCCC
CTGCTGAGGAAGAGACCAACCCCAAGGGCAGCGGCTGGTTGTACCATTGATGCCATCCGCACCTAC
CTGAACCTCATGGCAAGAGCAAGAAAGATGCTACCCTGGAGGCCTGTGCTGGTGCCTGCAGAACCTG
ACAGCCAGCAAGGGGCTGATGTCCAGTGGCATGAGCCAGTTGATTGGGCTGAAGGAAAAGGGCCTGCCA
CAAATTGCCCGCCTCCTGCAATCTGGCAACTCTGATGTGGTGGGTCGGTCCGGAGCCTCCCTCCTGAGAAC
ATGTCCCGCCACCTCTGCTGCACAGAGTGTGGGGAACCAAGGTGTTCCCGGAGGTGACCAGGCTCCTC
ACCAGCCACACTGGCAATACCAGCAACTCCGAAGACATCTTGTCTCGGCCTGCTACACTGTGAGGAAC
CTGATGGCCTCGCAGCCACAACCTGGCAAGCAGTACTTCTCCAGCAGCATGCTCAACAACATCATCAAC
CTGTGCCGAAGCAGTGCCTCACCAAGGCCGAGAAGTGGCCGGCTTCTCTGTCTGACATGTGGTCC
AGCAAGGAACTGCAGGGTGTCTCAGACAGCAAGGTTTCGATAGGAACATGCTGGAACTTAGCTGGG
GCCAACAGCCTCAGGAACTTACCTCCCGATTCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_000299

Insert Size: 2244 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_000299.2
RefSeq Size:	5447 bp
RefSeq ORF:	2244 bp
Locus ID:	5317
UniProt ID:	Q13835
Cytogenetics:	1q32.1
Protein Families:	Druggable Genome
MW:	82.9 kDa
Gene Summary:	<p>This gene encodes a member of the arm-repeat (armadillo) and plakophilin gene families. Plakophilin proteins contain numerous armadillo repeats, localize to cell desmosomes and nuclei, and participate in linking cadherins to intermediate filaments in the cytoskeleton. This protein may be involved in molecular recruitment and stabilization during desmosome formation. Mutations in this gene have been associated with the ectodermal dysplasia/skin fragility syndrome. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2010]</p> <p>Transcript Variant: This variant (1b) represents the longer transcript variant and encodes the longer protein (isoform 1b). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>