

## Product datasheet for **SC300067**

### Cytokeratin 2 (KRT2) (NM\_000423) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Cytokeratin 2 (KRT2) (NM_000423) Human Untagged Clone
Tag:	Tag Free
Symbol:	KRT2
Synonyms:	CK-2e; K2e; KRT2A; KRT2E; KRTE
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC300067 representing NM\_000423.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGAGTTGTCAGATCTCTTGCAAATCTCGAGGAAGAGGAGGAGGTGGAGGAGGATTCGGGGCTTCAGC
AGCGGCTCAGCTGTGGTGTCTGGTGAAGCCGGAGATCAACTCCAGCTTCTCTGCTTGAGCCGCCAT
GGTGGTGGTGGCGGGGGCTTCGGTGGAGGCGGCTTTGGCAGTCGGAGTCTTGTGGCCTTGAGGGACC
AAGAGCATCTCCATTAGTGTGGCTGGAGGAGGTGGTGGCTTTGGCGCCGCTGGTGGATTTGGTGGCAGA
GGAGGTGGTTTTGGAGGCGGCAGCAGCTTTGGAGGTGGCAGCGGCTTCAGTGGTGGTGGTTTTGGTGGG
GGCGGCTTTGGTGGAGGCGGCTTTGGAGTTTTGGGGCCCTGGTGGTGGTGGAGTTTTAGGGGTCTT
GGTGGCTTTGGGCTGGAGGATACCCTGGTGGCATCCACGAAGTCTCTGTCAACCAGAGCCTCTGCAG
CCTCTCAACGTGAAAGTTGACCCAGAGATCCAGAATGTGAAGGCCAAGAGCGTGAGCAGATCAAACCT
CTCAACAACAAATTTGCCTCCTTATTGACAAGGTGCGGTTCTTGGAGCAGCAGAACCAGGTGTTACAG
ACCAAATGGGAGCTGCTACAACAATGAATGTTGGCACCCGCCCATCAACTGGAGCCCATCTTCCAG
GGGTATATCGACAGCCTCAAGAGATATCTGGATGGGCTCACTGCAGAAAGAACATCACAGAAATTCAGAG
CTGAATAACATGCAGGATCTTGTGGAGGATTATAAGAAGAAGTATGAGGATGAAATCAATAAGCGCACA
GCTGCTGAGAATGATTTTGTGACGCTTAAAAAGGACGTGGACAATGCCTACATGATAAAGGTGGAGTTG
CAGTCCAAGGTGGACCTGCTGAACCAGGAAATTGAGTTTTCTGAAAGTTCTCTATGATGCGGAGATATCC
CAGATACATCAGAGTGTCACTGACACCAACGTCATCCTCTCCATGGACAACAGCCGAACCTGGACTTG
GATAGCATCATCGCCGAGGTCAAGGCCAGTATGAGGAGATCGCCAGAGGAGCAAGGAAGAAGCGGAG
GCCCTGTACCACAGCAAGTATGAGGAGCTCCAGGTGACTGTGGGAGACATGGAGACAGCCTGAAAGAG
ATCAAGATAGAGATCAGCGAGCTGAACCGCGTGATCCAGAGGCTGCAGGGGGAGATCGCACATGTGAAG
AAGCAGTGTAAGAATGTGCAAGATGCCATCGCAGATGCCGAGCAGCGTGGGGAGCATGCCCTCAAGGAT
GCCAGGAACAAGTTGAATGACCTGGAGGAGCCCTGCAGCAGGCCAAGGAGGACTTGGCGCGGCTGCTG
CGTGACTACCAGGAGCTGATGAACGTGAAGCTGGCCCTAGATGTGGAGATCGCCACCTACCGCAAACCTG
CTGGAGGGCGAGGAGTGCAGGATGTCTGGAGACCTCAGCAGCAATGTGACTGTGTCTGTGACAAGCAGC
ACCATTTTCATCAAATGTGGCATCCAAGGCTGCCTTTGGAGGTTCTGGAGGTAGAGGGTCCAGTTCGGGA
GGAGGATACAGCTCTGGAAGCAGCAGTTATGGCTCTGGAGGCCGACAGTCTGGCTCCAGAGGCGGTAGT
GGAGGAGGAGTTCTATCTCTGGAGGAGGATATGGCTCTGGCGGTGGTTCTGGAGGAAGATACGGATCT
GGTGGTGGCTCTAAGGGAGGGTCCATCTCTGGAGGAGGATATGGCTCTGGAGGTGAAAAACACAGCTCT
GGAGGTGGCTCTAGAGGAGGCTCCAGCTCTGGAGGAGGATATGGCTCTGGAGGTGGGGGTTCTAGCTCT
GTAAGGGTAGCTCAGGTGAAGCTTTTGGTTCAGCGTGACCTTCTCTTTTAGATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_000423

**Insert Size:** 1920 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_000423.2](#)

RefSeq Size: 2459 bp

RefSeq ORF: 1920 bp

Locus ID: 3849

UniProt ID: [P35908](#)

Cytogenetics: 12q13.13

MW: 65.4 kDa

**Gene Summary:** The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This type II cytokeratin is expressed largely in the upper spinous layer of epidermal keratinocytes and mutations in this gene have been associated with bullous congenital ichthyosiform erythroderma. The type II cytokeratins are clustered in a region of chromosome 12q12-q13. [provided by RefSeq, Jul 2008]