

Product datasheet for **SC300076**

AMPD3 (NM_000480) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AMPD3 (NM_000480) Human Untagged Clone
Tag:	Tag Free
Symbol:	AMPD3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC300076 representing NM_000480.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCCCTGTCGTCGGAACCCGCTGAGATGCCCGGCCAGTTTCCCAAGCTGAACATCTCTGAAGTGGAT
GAGCAAGTCCGGCTCTGGCGGAGAAGGTGTTTGTAAAGTGTCCGAGAAGAGGACAGCAAAGATGCC
CTGTCCCTGTTCACTGTCCAGAGGACTGCCCATCGGGCAAAAGGAAGCCAAGGAGAGGGAGCTGCAG
AAGGAGCTGGCAGAGCAGAAGTCTGTGGAGACCGCAAAAAGAAAAGTTTCAAGATGATTTCGGTCC
CAGTCCCTGTCTCTGCAATGCCGCCACAGCAAGATTGGAAGGGCCCCCGGCAGCCAGTCCGGCCATG
TCTCCACAACCCCTGTGGTCACTGGAGCCACTTCCCTGCCACGCCAGCACCCCTATGCCATGCCTGAG
TTCCAGCGGGTACCATCAGCGGAGATTACTGTCCGGGATCACTTTGGAGGACTATGAGCAGGCAGCC
AAGAGTCTGGCCAAGGCCCTAATGATCCGGGAGAAGTATGCGCGGCTCGCCTACCACCGCTTCCCGCGG
ATCACATCCAGTACCTGGTTCATCCGCGGGCGGATACTGCACCTCCGGAAGAGGGCCTTCCAGACTTC
CACCTCTCCACTGCCCAAGAACCCCTACTGCCTGGATGATGCACCCCCAACCTGGATTACTTG
GTCACATGCAGGGGGCATCCTTTTGTGTATGATAACAAGAAGATGTGGAGCACCAGGAGCCGCAC
AGCCTACCCTACCCGACCTGGAGACCTACACGGTGGACATGAGCCACATCCTGGCTCTATCACCGAT
GGCCCCACGAAAACCTATTGTCACCGGCGACTGAACTTTCTGGAATCCAAGTTGAGCCTTCATGAGATG
TTAAACGAAATGTCCGAGTTCAAAGAGTTGAAGAGTAACCCCAACCGGGACTTCTATAACGTGAGAAAG
GTGGACACACACATCCATGCGGCCGCTGCATGAACAAAAGCATCTGCTGCGTTCATCAAGCACACA
TACCAGACGGAGCCTGACAGGACTGTGGCAGAGAAGCGGGGCCGGAAGATCACCTGCGGCAGGTGTTT
GACCGCTGCACATGGACCCCTACGACCTCACTGTGGACTCACTGGATGTCCACGCGGGCCGCGCAGACA
TCCACCGCTTGGACAAGTTCACTCAAATACAACCTGTGGGGCCAGTGGAGTGCCTGACCTGTAT
TTGAAAACCTGAAAACCTATCTGGGAGGAGAGTACTTTGCTCGGATGGTCAAGGAGGTTGCCCGGAGCTG
GAGGAGAGCAAGTACCAGTACTCAGAGCCACGGCTCTCCATCTACGGCCGAGTCTGAGGAGTGGCCC
AACCTGGCTACTGGTTCATCCAGCACAAAGTCTACTCTCCAACATGCGCTGGATCATCCAGGTGCC
CGGATTTATGACATATTTAGGTCAAAGAAGCTGCTGCCAACTTTGGGAAGATGCTGGAGAACATCTTC
CTGCCCTTTTCAAGGCCACTATCAACCCCAAGATCATCGAGAGCTTACCTCTTCTTAAATATGTG
ACGGGGTTTGACAGCGTGGATGATGAGTCCAAGCACAGCGACCACATGTTTTCCGACAAGAGCCCAAAC
CCGGACGTCTGGACCAGTGAAGCAGAACCCACCCTACAGCTACTACCTGTACTACATGATGCCAACATC
ATGGTGTCAACAACCTCCGACGGAGCGGGCCTGAGCAGTTCCTGTTCCGGCCGACTGTGGGGAA
GCCGGCTCCATCACCCACCTGGTGTCTGCTTCTCACTGCTGACAACATTTCCACGGGCTGCTCCTC
AAGAAGAGTCCGGTATTGCAATCTCTACTACCTTGCTCAGATCCCCATTGCCATGTCTCCTTATGC
AACACAGTTTGTTCCTCGAATATTTCAAAGAACCCTCTGAGGGAATTCCTACACAAGGGACTGCATGTT
TCTCTTTCCACCGATGACCCCATGCAGTTCCTACTACAGGAAGGAAGCACTTATGGAAGAATAGCCATT
GCAGCTCAAGTGTGGAAGCTGAGCACCTGCGACCTGTGTGAGATCGCCAGGAACAGCGTGTGCAGAGC
GGCCTCTCGCATCAGGAAAAGCAAAAGTTTCTGGGACAAAATTATTATAAAGAAGGACCTGAAGGAAAT
GATATTCGAAAAGCAAAATGTGGCTCAGATCCGGATGGCATTCCGATATGAGACCTTATGCAATGAGCTC
AGTTCCTGTCTGATGCTATGAAATCAGAAGAGATCACCGCCTTGACCAACTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_000480

Insert Size: 2331 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_000480.2](#)

RefSeq Size: 4371 bp

RefSeq ORF: 2331 bp

Locus ID: 272

UniProt ID: [Q01432](#)

Cytogenetics: 11p15.4

Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Purine metabolism
MW:	89.7 kDa
Gene Summary:	<p>This gene encodes a member of the AMP deaminase gene family. The encoded protein is a highly regulated enzyme that catalyzes the hydrolytic deamination of adenosine monophosphate to inosine monophosphate, a branch point in the adenylate catabolic pathway. This gene encodes the erythrocyte (E) isoforms, whereas other family members encode isoforms that predominate in muscle (M) and liver (L) cells. Mutations in this gene lead to the clinically asymptomatic, autosomal recessive condition erythrocyte AMP deaminase deficiency. Alternatively spliced transcript variants encoding different isoforms of this gene have been described. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (1A).</p>